


**INNOVATIVE**  
**Medical Solutions**  
Specialists in Medical Electronics

# IBP Interface Cable and Transducers

1 July 2015

Phone: 1300 884 612  
Fax: 1300 884 613  
Email: [sales@innovativemedical.com.au](mailto:sales@innovativemedical.com.au)



# 1 IBP Adapter cable for Transducer

## 有创血压电缆

Compatible with major patient monitor brands

兼容主要品牌

- |                   |        |
|-------------------|--------|
| Datex-Ohmeda      | 德恩-欧美达 |
| Datascope-Mindray | 迈瑞     |
| Marquette         | 马奎     |
| Nihon kohden      | 日本光电   |
| Philips           | 飞利浦    |
| Spacelabs         | 太空     |
| Siemens           | 西门子    |
| ...               | ....   |

Specification  
产品规格

Cable material 电缆物料	TPU Jacket, Gray TUP料, 灰色
Connector 连接器	Nylon plug & Golden plated pins 尼龙插座和镀金针
Cable length 电缆长度	4.0m/0.3m
Cable Diameter 电缆直径	5.00mm

Available cables for disposable transducer  
用于一次性传感器的电缆

- Argon (AG)  
阿尔贡
- B.Braun (BB)  
贝朗
- B.D. (BD)  
贝达
- Edwards (ED)  
爱德华
- Medex/Abbort (MX)  
美代/雅培
- Utah (UT)  
尤它
- PVB



IBP Adapter cables

IBP Trans ducer Adapter

IBP Trans ducer

Reference Number Rule Description  
产品型号说明

BC-HP-UT/2	
BC	IBP cable 有创压电缆
HP	Factory/model (see others above) 厂家/型号 (其他, 看如上所示)
UT	Utah( transducer type) 尤它 (传感器类型),其他型号见如上示
2	the number of channel, default is single channel 通道数; 不写通道数, 则默认为单通道

Available cables for reusable transducer  
用于可重复使用传感器的电缆

Medex Logical/MX960 (美代)



Transducer Type 传感器类型	Argon 阿尔贡	B.Braun 贝朗	BD 贝达	Edward 爱德华	Medex / Abbott 美代/雅培	Utah 尤它	Medex Logical 美代	PVB	Cable Length 电缆长度	Cable Connector on Monitor end 仪器端电缆接头
<b>Cable Connector on transducer end 传感器端电缆接头</b>										
<b>Criticare</b>										
Model 2200	BC-CSI-AG C3	BC-CSI-BB C3	BC-CSI-BD C8	BC-CSI-ED C3	BC-CSI-MX C3	BC-CSI-UT C3	BC-CSI-MX1 D1	BC-CSI-PVB C3	4.0m	
<b>Colin</b>										
	BC-CL-AG C3	BC-CL-BB C3	BC-CL-BD C8	BC-CL-ED C3	BC-CL-MX C3	BC-CL-UT C3	BC-CL-MX1 D1	BC-CL-PVB C3	4.0m	
<b>Biolight</b>										
V6 Vital Sign View A2/A3/A5/A6	BC-BLT-AG C3	BC-BLT-BB C3	BC-BLT-BD C8	BC-BLT-ED C3	BC-BLT-MX C3	BC-BLT-UT C3	BC-BLT-MX1 D1	BC-BLT-PVB C3	4.0m	
<b>GE-Datex</b>										
Airshield, Burdick, Kone	BC-DX-AG C3	BC-DX-BB C3	BC-DX-BD C8	BC-DX-ED C3	BC-DX-MX C3	BC-DX-UT C3	BC-DX-MX1 D1	BC-DX-PVB C3	4.0m	
<b>GE-Hellige</b>										
	BC-HL-AG C3	BC-HL-BB C3	BC-HL-BD C8	BC-HL-ED C3	BC-HL-MX C3	BC-HL-UT C3	BC-HL-MX1 D1	BC-HL-PVB C3	4.0m	
<b>GE-Marquette</b>										
Eagle/dash1000/2000 3000/3100/4000/Solar 7000/8000/9000/Tram Module 100-600, Cath-Lab, MAC-Lab Dash 1/Tram(all versions) E-PSMP/E-PRESTN/PDM	BC-MQ-AG C3	BC-MQ-BB 22339301 C3	BC-MQ-BD 2016995-001 C8	BC-MQ-ED 2021197-001 C3	BC-MQ-MX 2021196-001 C3	BC-MQ-UT 2016997-001 C3	BC-MQ-MX1 D1	BC-MQ-PVB C3	4.0m	
	BC-MQ-AG/2 C8	BC-MQ-BB/2 C8	BC-MQ-BD/2 2016995-003 D3	BC-MQ-ED/2 2021197-003 C8	BC-MQ-MX/2 2021196-003 C8	BC-MQ-UT/2 2016997-003 C8	BC-MQ-MX1/2 D4	BC-MQ-PVB/2 C8		
<b>Mostcare</b>										
	BC-MC-AG C3	BC-MC-BB C3	BC-MC-BD C8	BC-MC-ED C3	BC-MC-MX C3	BC-MC-UT C3	BC-MC-MX1 D1	BC-MC-PVB C3	4.0m	
<b>Mennen</b>										
Envoy, vitalogik 5500	BC-MN-AG C3	BC-MN-BB C3	BC-MN-BD C8	BC-MN-ED C3	BC-MN-MX C3	BC-MN-UT C3	BC-MN-MX1 D1	BC-MN-PVB C3	4.0m	
<b>Mindray/Datascope</b>										
Beneview T5/T6/T8 DPM 6, DPM 7	BC-MR-AG C3	BC-MR-BB C3	BC-MR-BD 040-000053-00 C8	BC-MR-ED 040-000054-00 C3	BC-MR-MX 040-000052-00 C3	BC-MR-UT C3	BC-MR-MX1 D1	BC-MR-PVB C3	4.0m	
<b>Datascope</b>										
Passport V, V12, V21	BC-DT1-AG C3	BC-DT1-BB C3	BC-DT1-BD 001C-30-70757 C8	BC-DT1-ED 040-000053-00 C3	BC-DT1-MX 040-000054-00 C3	BC-DT1-UT 040-000052-00 C3	BC-DT1-MX1 D1	BC-DT1-PVB C3	4.0m	
<b>Philips/HP</b>										
78352A, M1165A, Meriv/Vrds 24, CMS 24,M3/M4,M1175A, M1350B,M1205A, M1350A,M1353A, 8030A,8040A	BC-HP-AG C3	BC-HP-BB M1634A C3	BC-HP-BD C8	BC-HP-ED C3	BC-HP-MX C3	BC-HP-UT C3	BC-HP-MX1 D1	BC-HP-PVB C3	4.0m	
<b>Siemens</b>										
Sirecut 400-1281, SC8000	BC-SM1-AG C3	BC-SM1-BB C3	BC-SM1-BD C8	BC-SM1-ED C3	BC-SM1-MX C3	BC-SM1-UT C3	BC-SM1-MX1 D1	BC-SM1-PVB C3	4.0m	
	BC-SM1-AG/2 C8	BC-SM1-BB/2 C8	BC-SM1-BD/2 D3	BC-SM1-ED/2 C8	BC-SM1-MX/2 C8	BC-SM1-UT/2 C8	BC-SM1-MX1/2 D4	BC-SM1-PVB/2 C8		
SC 6000,7000,9000	BC-SM2-AG C3	BC-SM2-BB C3	BC-SM2-BD C8	BC-SM2-ED C3	BC-SM2-MX C3	BC-SM2-UT C3	BC-SM2-MX1 D1	BC-SM2-PVB C3	4.0m	
	BC-SM2-AG/2 C8	BC-SM2-BB/2 C8	BC-SM2-BD/2 D3	BC-SM2-ED/2 C8	BC-SM2-MX/2 C8	BC-SM2-UT/2 C8	BC-SM2-MX1/2 D4	BC-SM2-PVB/2 C8		

Transducer Type 传感器类型	Argon 阿尔贡	B.Braun 贝朗	BD 贝达	Edward 爱德华	Medex / Abbott 美代/雅培	Utah 尤它	Medex Logical 美代	PVB	Cable Length 电缆长度	Cable Connector on Monitor end 仪器端电缆接头
<b>Cable Connector on transducer end 传感器端电缆接头</b>										
<b>Drager</b>										
Cato, Cicero EM, Vitara Parameterbox PB 8800; PM8010; PM8014; PM 8040; PM 8060, Sulla, UM3, UM3.1	BC-DG-AG C3	BC-DG-BB C3	BC-DG-BD C8	BC-DG-ED C3	BC-DG-MX C3	BC-DG-UT C3	BC-DG-MX1 D1	BC-DG-PVB C3	4.0m	
Infinity, Gamma X Xi Delta/delta Xi Kappa/kappa Xi Vista Xi	BC-S-AG/2 D5	BC-S-BB/2 D5	BC-S-BD/2 D5	BC-S-ED/2 D5	BC-S-MX/2 D5	BC-S-UT/2 D5	BC-S-MX1/2 D6	BC-S-PVB/2 D5	4.0m	
<b>Datascope</b>										
System 80/82/83/84 90/90T/95/98/97IABP 850/870/P2/P3/2000 2000A/2001/2001A 2002/2002A/22001 3000/Passport	BC-DT-AG C3	BC-DT-BB C3	BC-DT-BD 0012-00-1245 C8	BC-DT-ED C3	BC-DT-MX C3	BC-DT-UT C3	BC-DT-MX1 D1	BC-DT-PVB C3	4.0m	
<b>Fukuda Denshi</b>										
DYNASCOPE/Expert DS5100E/DS5300W DS5400	BC-FD-AG C3	BC-FD-BB C3	BC-FD-BD C8	BC-FD-ED C3	BC-FD-MX C3	BC-FD-UT C3	BC-FD-MX1 D1	BC-FD-PVB C3	4.0m	
<b>Nihon Kohden</b>										
AP-311P -800PA BMS 2100-8500A0EC 3200-8120A/OMP-3221E ~n7203/Mercury 9510K	BC-NK1-AG C3	BC-NK1-BB C3	BC-NK1-BD C8	BC-NK1-ED C3	BC-NK1-MX C3	BC-NK1-UT C3	BC-NK1-MX1 D1	BC-NK1-PVB C3	4.0m	
New	BC-NK2-AG C8	BC-NK2-BB C8	BC-NK2-BD JP-900P D3	BC-NK2-ED JP-920P/3.5m JP-921P/1.5m C8	BC-NK2-MX C8	BC-NK2-UT JP-960P C8	BC-NK2-MX1 D1	BC-NK2-PVB C8	4.0m	
<b>Kontron</b>										
Minimon 6/9/12/14 Fetalmon 2A 2BFetalogic 4000 Fetasonde 5/5A/5B K2000 Ballon Pump KAAT II/Module 10-7286 Supernon	BC-KTN-AG C3	BC-KTN-BB C3	BC-KTN-BD C8	BC-KTN-ED C3	BC-KTN-MX C3	BC-KTN-UT C3	BC-KTN-MX1 D1	BC-KTN-PVB C3	4.0m	
<b>Spacelabs</b>										
Ultraview 90496, 91496	BC-SL-AG C5	BC-SL-BB C5	BC-SL-BD D1	BC-SL-ED 700-0293-00 C5	BC-SL-MX 700-0295-00 C5	BC-SL-UT C5	BC-SL-MX1 D1	BC-SL-PVB C5	4.0m	
	BC-SL-AG/2 D1	BC-SL-BB/2 D1	BC-SL-BD/2 D3	BC-SL-ED/2 700-0292-00 D1	BC-SL-MX/2 700-0294-00 D1	BC-SL-UT/2 D1	BC-SL-MX1/2 D4	BC-SL-PVB/2 D1		
<b>Spacelabs &amp; General</b>										
BCI, IVY,Citikon, Invivo,Mennen,NDE, PB,Nellcor,Tektronix, Viatak,Mindray	BC-6P-AG B8	BC-6P-BB B8	BC-6P-BD C8	BC-6P-ED B8	BC-6P-MX B8	BC-6P-UT B8	BC-6P-MX1 D1	BC-6P-PVB B8	4.0m	
<b>Stockert</b>										
	BC-SK-AG C8	BC-SK-BB C8	BC-SK-BD D3	BC-SK-ED C8	BC-SK-MX C8	BC-6P-UT C8	BC-SK-MX1 D1	BC-SK-PVB C8	4.0m	

# ② IBP Convert Cable

## 有创电缆转接线

### Reference Number Rule Description 产品型号说明

BC-SL-AD/2	
BC	IBP cable 有创血压电缆
SL	Factory/model (see others on previous page) 厂家/型号 (其他参考前页)
AD	Convert Cable.AD: Converted to AAMI Adapter(MQ: to marquette;SM1 or SM2: to Siemens) 转接线 AD: 连接6针通用有创电缆 (MQ: 连接marquette电缆; SM1或SM2: 连接西门子电缆)
2	2-channel,default is single channel 2-通道数; 不写通道数, 则默认为单通道



Cable end Connector 电缆端接头	Unimed REF# 尤迈型号	OEM PN 原厂型号	Used with 连接电缆	Reference 备注	Cable Length 电缆长度	Monitor End Connector 仪器端接头	Price code 价格代码
<b>GE Marquette</b>							
	BC-MQ-AD/2		6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous for Marquette 参考前页面展示的马赛型号	0.3m		D3
	BC-MQ-MQ/2	2005772-001	Marquette IBP adapter Cable 马赛有创血压转接线	See the previous for Marquette 参考前页面展示的马赛型号	0.3m		D5
<b>Spacelabs</b>							
	BC-SL-AD	700-0791-00	6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous page for Spacelabs 参考前页面展示的太空型号	0.3m		D1
	BC-SL-AD/2	700-0028-00	6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous page for Spacelabs 参考前页面展示的太空型号	0.3m		D5
<b>Datascopes/Mindray</b>							
	BC-MR-AD	0010-20-42795	6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous page for Mindray. 参考前页面展示的心迈型号	0.3m		C3
	BC-DT1-AD	040-000096-00	6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous page for Mindray 参考前页面展示的心迈型号	0.3m		C3
	BC-DT1-DT		6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous page for Mindray & Datascopes 参考前页面展示的心迈 & 达泰型号	0.3m		C3
<b>Siemens/Drager</b>							
	BC-S-SM1/2	5731281	SM1 conn. 2XIBP Cable SM1接头 2根有创血压电缆	See the previous Page for Siemens 参考前页面展示的西子型号	0.3m		D6
	BC-S-SM2/2	5592147	SM2 conn. 2XIBP Cable SM2接头 2根有创血压电缆	See the previous Page for Siemens 参考前页面展示的西子型号	0.3m		D6
	BC-SM2-SM1/2	5588095 Drager Gamma XL and Delta XL	SM1 conn. 2XIBP Cable SM1接头 2根有创血压电缆	See the previous Page for Siemens 参考前页面展示的西子型号	0.3m		D5
	BC-SM1-SM2		SM2 conn. IBP Cable 有创血压电缆	See the previous Page for Siemens 参考前页面展示的西子型号	0.3m		C3
	BC-SM2-SM1	3368383	SM1 conn. IBP Cable 有创血压电缆	See the previous Page for Siemens 参考前页面展示的西子型号	0.3m		C3

### ③ IBP Transducer Adapter 有创血压传感器转接线




Reference Number Rule Description  
产品型号说明

BC-MX-AG	
BC	IBP cable有创血压电缆
MX	Factory/model, used with Medex adaptor Cable厂家/型号, 与美代有创血压电缆配套。
AG	Factory/model, used with Argon disposable transducer厂家/型号, 与阿尔贡一次性传感器配套。 (Factory/model see others above)(具体厂家/型号, 参考前面页的描述)

Used with brands IBP Adapter cable & Major types disposable transducer  
连接各品牌IBP电缆及主要一次性有创压传感器

Transducer Type 传感器类型	Argon 阿尔贡	B.Braun 贝朗	BD 贝达	Edward 爱德华	Medex / Abbott 美代/雅培	Utah 尤它	Medex Logical 美代	Cable Length 电缆长度
With transducer conn. 与传感器接头	NA		NA	NA			NA	
With IBP cable conn. 与电缆接头								
		BC-BB-AG C3	BC-BD-AG C5	BC-ED-AG C3	BC-MX-AG C3	BC-UT-AG C3	BC-ML-AG C3	0.3m
	BC-AG-BB C3		BC-BD-BB C5	BC-ED-BB C3	BC-MX-BB C3	BC-UT-BB C3	BC-ML-BB C3	0.3m
	BC-AG-BD C3	BC-BB-BD C3		BC-ED-BD C3	BC-MX-BD C3	BC-UT-BD C3	BC-ML-BD C3	0.3m
	BC-AG-ED C3	BC-BB-ED C3	BC-BD-ED C5		BC-MX-ED C3	BC-UT-ED C3	BC-ML-ED C3	0.3m
	BC-AG-MX C3	BC-BB-MX C3	BC-BD-MX C5	BC-ED-MX C3		BC-UT-MX (97005-C) C3	BC-ML-MX C3	0.3m
	BC-AG-UT C3	BC-BB-UT C3	BC-BD-UT C5	BC-ED-UT C3	BC-MX-UT C3		BC-ML-UT C3	0.3m
	BC-AG-MX1 C3	BC-BB-MX1 C3	BC-BD-MX1 C8	BC-ED-MX1 C3	BC-MX-MX1 C3	BC-UT-MX1 C3		0.3m



### ④ Disposable Transducer 一次性有创压传感器

