



**INNOVATIVE**  
**Medical Solutions**  
Specialists in Medical Electronics

# Smart Temp Probes

1 July 2015

Phone: 1300 884 612

Fax: 1300 884 613

Email: [sales@innovativemedical.com.au](mailto:sales@innovativemedical.com.au)



### 3.1 Smart Temp. Probe 快速体温探头



Reference Number Rule Description  
产品型号说明

TOP-AR-001
T Temperature 体温
O Oral口腔 (Rectal直肠)
P Probe 探头
AR Factory Code 厂家代码 ( AR=Alaris, GE=GE Healthcare, MR=Mindray)
001 Factory /model 厂家/型号

### ③ Smart Temp. Probe 快速体温探头

Compatible with major patient monitor brands 兼容主要品牌

Major brands temperature monitoring system  
主要品牌体温监护系统

#### Features 产品特征

- Accurate measurement 精确测量
- 15s Quick response 15秒快速读取数据
- Low noise signal transfer 抗干扰性强
- Fully compatible with OEM technology 完全兼容原厂技术
- 12 month warranty 12个月保修期

#### Specification 产品规格

Probe type 探头类型	Oral;Rectal 口腔; 直肠
Cable length 电缆长度	4m unless requirement applied 4米 特殊要求除外
Resistor type 电阻类别	NTC NTC
Cable Material 线材材质	TPU jacket,Light gray TPU外皮, 灰色
Package 包装	1pcs per bag 每袋一包

OEM & Models 原厂品牌 / 机型	Unimed REF# 尤迈型号	OEM REF# OEM型号	Application 应用	Cable length 线长(米)	Monitor end plug 机器端插头	Price code 价格代码
Alaris						
VITAL CHECK 4400/4410/4500	TOP-AR-001	3880L	Oral 成人口腔	4m		D5
	TRP-AR-001	3882L	Rectal成人直肠	4m		D5
Turbo Temp 2180C/2180CX01EE	TOP-AR-002	2885C	Oral 成人口腔	4m		D5
	TRP-AR-002	2886C	Rectal成人直肠	4m		D5
Turbo Temp 2185B/2185BX01EE	TOP-AR-003	2887B	Oral 成人口腔	4m		D5
	TRP-AR-003	2888B	Rectal成人直肠	4m		D5
GE Healthcare						
DINAMAP Pro200/400 V1	TOP-AR-001	088012	Oral 成人口腔	4m		D5
	TRP-AR-001	088016	Rectal成人直肠	4m		D5
DINAMAP Pro200/400 V2 ProCare200/400	TOP-AR-002	2008774-001	Oral 成人口腔	4m		D5
	TRP-AR-002	2008775-001	Rectal成人直肠	4m		D5
DINAMAP 9300 Series DINAMAP XL	TOP-GE-001	008975,008978	Oral 成人口腔	4m		D5
	TRP-GE-001	008976	Rectal成人直肠	4m		D5
Mindray/迈瑞						
Accutorr Plus	TOP-MR-001	0206-00-0725-01	Oral 成人口腔	4m		D5
	TRP-MR-001	0206-00-0725-02	Rectal成人直肠	4m		D5