

INNOVATIVE
Medical Solutions
Specialists in Medical Electronics

Reusable SpO₂ Sensors


1 July 2015

Phone: 1300 884 612
Fax: 1300 884 613
Email: sales@innovativemedical.com.au



Compatible with major patient monitor brands

- Nellcor
- Ohmeda
- Philips
- Marquette
- Spacelabs
- Novamatrix
- Criticare
- BCI
- Datascope
- Datex
- Nihon Kohden
- Palco



1 Reusable SpO₂ Sensor

Specification

Sensor Connector	OEM compatible, Nylon adapter with gold plated pins
LED & Photo Detector	660nm/905nm/940nm plastic or ceramic frame LED; 4mm or 8mm size plastic flame PD
Cable material & color	TPU jacket; light gray or bright blue
Cable length & diameter	1.1m/1.5m/3.0m; φ 3.20mm/4.00mm or 5.00x2.50mmtwo line cable
Sensor house	ABS shell(finger clip); Silicon rubber boots(soft and multi-site)

Series	U400	U100	U400S	U100S	U200S	U300
REF to Fig	A	B	C	D	E	F
Name	Adult Finger	Pediatric Finger	Adult soft tip	Pediatric Soft tip	Infant Soft tip	Neonate Wrap
Patient size	+40kg	15-40kg	+40kg	15-40kg	3-15kg	1-3kg
Cable color	Gray Cable	Gray Cable	Blue Cable	Gray two line	Gray two line	Gray two line

Series	U800	U900	U600	U600B	U800E	U800V
REF to Fig	G	H	I	J	K	L
Name	Multi-site	Adult Ear clip	Adult/Neo. Wrap	Pedia/Inf. Wrap	Multi-site	Forehead/Rectal
Patient size	All patients	+40kg	+40kg or -3kg	15-40kg	All patients	Adult
Cable color	Gray two line		Gray two line	Gray two line	Sensormate Used with DPT tapes	Gray Cable

How to identify the REF# of Unimed Spo2 sensors

Reference Number Rule Description

U403S-01XE	
U	Unimed brand
4	adult sensor (1: pediatric; 2: infant; 3: neonate; 5: disposable; 6 wrap; 7: ext-cable; 8: multi-site; 9: ear clip)
03	1.1m cable length (03: 1.1m; 05: 1.5m; 08: 2.2m; 10: 3.0m)
S	Soft sensor (B:U style shell; Default: shell sensor or neonate wrap, alternative wrap, multi-site)
01	Nellcor code (02:Ohmeda; 03:Novamatrix; 04: Spacelabs; 05: Criticare; 06: BCI...)
X	Oximax technology (M: M technology, D: digital technology; Default: standard)
E	Sensormate sensor (Default: non-Sensormate sensor)

Cross Reference

Ending "X" can be replaced by "XP" or "P" for fully compatible with all oximax oximetry

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
3F Medical (China)								
IRIS 300/IRIS 400 Phoebe PMS8300/8310 PhoebeProPMS8300/ 8310PMS8110A			U410-189	Adult finger /3.0m	C7	No need		
			U410S-189	Adult soft /3.0m	C8			
			U810-189	Multi-site Y 3.0m	C7			
Adecon (China)								
MCF-21L, DK-8000D,DK-8000L			U410-120	Adult finger /3.0m	C7	No need		
			U410S-120	Adult soft /3.0m	C8			
			U810-120	Multi-site Y /3.0m	C7			
Anke (China)								
ASC-553			U410-37	Adult finger /3.0m	C7	No need		
			U410S-37	Adult soft /3.0m	C8			
			U810-37	Multi-site Y /3.0m	C7			
ASC-553 (*Oximax tech)			U410-78P	Adult finger /3.0m	C8	No need		
			U410S-78P	Adult soft /3.0m	D1			
			U810-78P	Multi-site Y /3.0m	C8			
ASC-553A			U410-82	Adult finger /3.0m	C7	No need		
			U410S-82	Adult soft /3.0m	C8			
			U810-82	Multi-site Y /3.0m	C7			
ASC-553A (*Ohmeda Tech)			U410-107	Adult finger /3.0m	C7	No need		
			U410S-107	Adult soft /3.0m	C8			
			U810-107	Multi-site Y /3.0m	C7			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.	
Artema/S&W									
Diascope,8000, Athena (*Non-Oximax)		DS100A	U403-01	Adult finger成人指夹/1.1m	C3	U708-80			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5				
	D-YS	U803-01	Multi-site Y型多功能/1.1m	C3					
			U410-80	Adult finger 成人指夹/3.0m	C8				No need 无需转接线
			U410S-80	Adult soft 成人软指套/3.0m	D1				
		U810-80	Multi-site Y型多功能/3.0m	C8					
S&W			U410-19	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-19	Adult soft 成人软指套/3.0m	D1				
			U810-19	Multi-site Y型多功能/3.0m	C8				
Baxter/Simed									
ASAT, S-50/100/100a			U410-12	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-12	Adult soft 成人软指套/3.0m	D1				
			U810-12	Multi-site Y型多功能/3.0m	C8				
Monarch/Explorer (*Oximax)		DS100A	U403-01	Adult finger成人指夹/1.1m	C3	U708-01			
		DS100A	U403S-01	Adult soft成人软指套/1.1m	C5				
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3				
BCI									
3300,3301,3302, 3303,3304,3401 (*BCI)		3044	U403-06	Adult finge 成人指夹/1.1m	C3	Depends U708-06 可选U708-06			
			U103-06	Pediatric finger 儿童指夹/1.1m	C5				
			U403S-06	Adult soft 成人软指套/1.1m	C5				
			U103S-06	Pediatric soft 儿童软指套/1.1m	C5				
			U203S-06	Infant soft 婴儿软指套/1.1m	C5				
			U303-06	Neonate wrap 新生儿捆绑/1.1m	C5				
	3043		U803-06	Multi-site Y型多功能/1.1m	C3				
			U603-06	Adult/Neo Wrap成人新生儿捆绑/1.1m	C3				
			U603B-06	Pedia/infant wrap/儿童婴儿捆绑/1.1m	C3				
			U903-06	Adult ear clip 成人耳夹/1.1m	C5				
3300,3301 6100,9100 Advisor, Autocorr, MiniTorr, Mini Torr plus (BCI*)		See above*	See the above sensor*	1.1m	See above sensors*	U708-76			
			U410-76	Adult finger 成人指夹/3.0m	C7				
			U410S-76	Adult soft 成人软指套/3.0m	C8				No need 无需转接线
			U810-76	Multi-site Y型多功能/3.0m	C7				

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Biosys								
BPM200 BPM700 Guardian		DS100A	U403-01	Adult finger /1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
			U410-27	Adult finger /3.0m	C8	No need		
			U410S-27	Adult soft /3.0m	D1			
			U810-27	Multi-site Y/3.0m	C8			
Biolight宝莱特(China)								
M6,M12,M69			U410-47	Adult finger /3.0m	C7	No need		
			U410S-47	Adult soft /3.0m	C8			
			U810-47	Multi-site Y/3.0m	C7			
M66,M69			U410-57	Adult finger /3.0m	C8	No need		
			U410S-57	Adult soft /3.0m	D1			
			U810-57	Multi-site Y/3.0m	C8			
(*Oximax)			U410-51P	Adult finger 3.0m	C8	No need		
			U410S-51P	Adult soft/3.0m	D1			
			U810-51P	Multi-site Y/3.0m	C8			
M700 (*Digital Tech)			U403-62D	Adult finger /1.1m	C8	No need		
			U403S-62D	Adult soft /1.1m	D1			
			U803-62D	Multi-site Y/1.1m	C8			
M8000,M9000,M8000A M9000A,BTD-352A, new M69 (*Digital Tech)			U410-62D	Adult finger /3.0m	D1	No need		
			U410S-62D	Adult soft /3.0m	D2			
			U810-62D	Multi-site Y/3.0m	D1			
M900,M7000,M9500 A8 Modular (*Digital Tech)			U410-68D	Adult finger /3.0m	D1	No need		
			U410S-68D	Adult soft /3.0m	D2			
			U810-68D	Multi-site Y/3.0m	D1			
Bionet								
		DS100A	U403-01	Adult finger /1.1m	C3	U708-58		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
BM3/BM5			U410-58	Adult finger /3.0m	C7	No need		
			U410S-58	Adult soft /3.0m	C8			
			U810-58	Multi-site Y/3.0m	C7			
BM3/BM5(New)			U410-108	Adult finger /3.0m	C7	No need		
			U410S-108	Adult soft /3.0m	C8			
			U810-108	Multi-site Y/3.0m	C7			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Biocare (China)								
BP-12B,BP-12C (*Digital Tech)			U410-61D	Adult finger /3.0m	D1	No need		
			U410S-61D	Adult soft /3.0m	D2			
			U810-61D	Multi-site Y/3.0m	D1			
BM-9000,BM-9000S FM-801,FM-801V (*Digital Tech)			U410-62D	Adult finger /3.0m	D1	No need		
			U410S-62D	Adult soft /3.0m	D2			
			U810-62D	Multi-site Y/3.0m	D1			
BM-9000			U410-29	Adult finger /3.0m	C7	No need		
			U410S-29	Adult soft/3.0m	C8			
			U810-29	Multi-site Y/3.0m	C7			
CAS Medical								
Oscillomate 9002 9303,9302S (*Nonin)		8000AA	U403-08	Adult finger /1.1m	C3	U708-01		
			U403S-08	Adult soft/1.1m	C5			
			U803-08	Multi-site Y/1.1m	C3			
740-750(M) (*M tech)		1863/ LNCS DCI	*U403-49	Adult finger /1.1m	C4	U708M-15		
		1863/ LNCS DCI	*U403S-49	Adult soft/1.1m	C6			
		2258/ LNCS Y1	*U803-49	Multi-site Y/1.1m	C4			
		LNOP DC-12	*U410-15	Adult finger /3.0m	C8	No need		
		LNOP DC-12	*U410S-15	Adult soft/3.0m	D1			
			*U810-15	Multi-site Y/3.0m	C8			
Charmcare								
CX100 (*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft/1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
ChoiceMed (China)								
MMED3001		DS100A	U403-01	Adult finger/1.1m	C3			
		DS100A	U403S-01	Adult soft/1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
MMED100H, MMED300A, MMED6000DP (*Digital tech)			U410-84D	Adult finger /3.0m	D1	No need		
			U410S-84D	Adult soft/3.0m	D2			
			U810-84D	Multi-site Y/3.0m	D1			
MMED6000DP (*BCI Tech)			U410-96	Adult finger /3.0m	C7	No need		
			U410S-96	Adult soft/3.0m	C8			
			U810-96	Multi-site Y/3.0m	C7			
Contec(China)								
			U410-29	Adult finger /3.0m	C7	No need		
			U410S-29	Adult soft 3.0m	C8			
			U810-29	Multi-site Y/3.0m	C7			
CMS60C/CW/D/DW,CMS60, PM50,PM60D,CMS5000, CMS5000C,CMS5100, CMS6000,CMS7000/ 8000/9000(*Digital Tech)			U410-46D	Adult finger /3.0m	C8	No need		
			U410S-46D	Adult soft/3.0m	D2			
			U810-46D	Multi-site Y/3.0m	D1			
Contec China								

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头						
Contec/康泰 (China)														
CMS50F/CMS50FW			U401-149	Adult finger 成人指夹/0.3m	C7	No need								
			U401S-149	Adult soft 成人软指套/0.3m	C8	无需转接线								
Colin/Omron科林 (欧姆龙)														
BP306, BP88, BP508 (Nellcor Non-Oximax) M570		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3									
			U403S-01	Adult soft 成人软指套/1.1m	C5	U708-01								
			D-YS	U803-01	Multi-site Y型多功能/1.1m	C3								
BP88S (*Non-Oximax)		See the above	See the above	See the above	See the above	U708-69								
									U410-69	Adult finger 成人指夹/3.0m	C8	No need		
									U410S-69	Adult soft 成人软指套/3.0m	D1	无需转接线		
BX 10NE (* Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5									
			U403S-01X	Adult soft 成人软指套/1.1m	C7	U710-70								
			D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5								
BX 10MA (* M tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4									
			*U403S-49	Adult soft 成人软指套/1.1m	C6	U708M-15								
			*U803-49	Multi-site Y型多功能/1.1m	C4									
			*U410-15	Adult finger 成人指夹/3.0m	C8	No need								
Corpuls		3044	U403-06	Adult finger 成人指夹/1.1m	C3									
			U403S-06	Adult soft 成人软指套/1.1m	C5	NA								
			U803-06	Multi-site Y型多功能/1.1m	C3									
08/16(M) (* Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3									
			U403S-01	Adult soft 成人软指套/1.1m	C5	NA								
			D-YS	U803-01	Multi-site Y型多功能/1.1m	C3								
08/16(M) M Tech			U410-152M	Adult finger 成人指夹/3.0m	C8	No need								
			U410S-152M	Adult soft 成人软指套/3.0m	D1	无需转接线								
			U810-152M	Multi-site Y型多功能/3.0m	C8									
Comen科曼 (China)														
Start-8000A/8000B/8000C Star-8000D/8000H Star-5000B/5000C, C50 (*Digital Tech)			U410-66D	Adult finger 成人指夹/3.0m	D1	No need								
			U410S-66D	Adult soft 成人软指套/3.0m	D2	无需转接线								
			U810-66D	Multi-site Y型多功能/3.0m	D1									
Start-8000E			U410-183	Adult finger 成人指夹/3.0m	C7	No need								
			U410S-183	Adult soft 成人软指套/3.0m	C8	无需转接线								
			U810-183	Multi-site Y型多功能/3.0m	C7									
C60, C80 (*Oximax Tech)			U410-187P	Adult finger 成人指夹/3.0m	D1	No need								
			U410S-187P	Adult soft 成人软指套/3.0m	D2	无需转接线								
			U810-187P	Multi-site Y型多功能/3.0m	D1									

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
Creative科瑞康 (China)								
UP-6000, PC-9000 (*Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3			
			U403S-01	Adult soft 成人软指套/1.1m	C5	U708-29		
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
			U410-29	Adult finger 成人指夹/3.0m	C7	No need		
			U410S-29	Adult soft 成人软指套/3.0m	C8	无需转接线		
UP-7000, UP-9000, PC-900, PC-900A, PC-900S, PC-900SN, PC-3000, PC-5000A/B, PC-9000B, PC-66A, PC-66B, PC-66V (*Digital Tech)			U810-29	Multi-site Y型多功能/3.0m	C7			
			U410-61D	Adult finger 成人指夹/3.0m	C8	No need		
			U410S-61D	Adult soft 成人软指套/3.0m	D1	无需转接线		
			U810-61D	Multi-site Y型多功能/3.0m	C8			
Criticare (CSI)								
510, 503, 504, 507, 508, 581, 1100, 8100, Poet, Scholar (CSI*)			934-10LN	U410-05	Adult finger 成人指夹/3.0m	C8		
			975AD-10	U110-05	Pediatric finger 儿童指夹/3.0m	D1		
				U410S-05	Adult soft 成人软指套/3.0m	D1	No need	
				U110S-05	Pediatric soft 儿童软指套/3.0m	D1	无需转接线	
			U210S-05	U310-05	Infant soft 婴儿软指套/3.0m	D1		
				U310-05	Neonate wrap 新生儿捆绑/3.0m	D1		
				U610-05	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8		
				U610B-05	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8		
				U810-05	Multi-site Y型多功能/3.0m	C8		
				U910-05	Adult ear clip 成人耳夹/3.0m	D1		
			U403-05	U403-05	Adult finger 成人指夹/1.1m	C3		
				U103-05	Pediatric finger 儿童指夹/1.1m	C4		
				U403S-05	Adult soft 成人软指套/1.1m	C4		
				U103S-05	Pediatric soft 儿童软指套/1.1m	C4		
				U203S-05	Infant soft 婴儿软指套/1.1m	C4	U708-05	
U303-75	Neonate wrap 新生儿捆绑/1.1m	C4						
U603-05	Multi-site Y型多功能/1.1m	C3						
U603B-05	Adult/Neo成人/新生儿捆绑/1.1m	C3						
U803-05	Pedia/infant wrap 儿童/婴儿捆绑/1.1m	C4						
U903-05	Adult ear clip 成人耳夹/1.1m	C4						
503SPOT, 504Plus, 504DX, 506DX (*CSI)			934-10DN	U410-75	Adult finger 成人指夹/3.0m	C7		
			U110-75	Pediatric finger 儿童指夹/3.0m	C8			
			U410S-75	Adult soft 成人软指套/3.0m	C8	No need		
			U110S-75	Pediatric soft 儿童软指套/3.0m	C8	无需转接线		
			U210S-75	Infant soft 婴儿软指套/3.0m	C8			
			U310-75	Neonate wrap 新生儿捆绑/3.0m	C8			
			U610-75	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C7			
			U610B-75	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C7			
			U810-75	Multi-site Y型多功能/3.0m	C7			
			U910-75	Adult ear clip 成人耳夹/3.0m	C8			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
Datascope								
Accusat, Accutor 3/4 sat Passport, Accutor plus (*Datascope)			U403-07	Adult finger 成人指夹/1.1m	C4	U708-07		
			U103-07	Pediatric finger 儿童指夹/1.1m	C6			
			U403S-07	Adult soft 成人软指套/1.1m	C6			
			U103S-07	Pediatric soft 儿童软指套/1.1m	C6			
			U203S-07	Infant soft 婴儿软指套/1.1m	C6			
			U303-07	Neonate wrap 新生儿捆绑/1.1m	C6			
			U603-07	Adult/Neo wrap 成人/新生儿捆绑/1.1m	C5			
			U603B-07	Pedia/infant wrap 儿童/婴儿捆绑/1.1m	C5			
			U803-07	Multi-site Y型多功能/1.1m	C5			
			U903-07	Adult ear clip 成人耳夹/1.1m	C4			
Accutor plus (*M Tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-15		
		1863/ LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6			
		2258/ LNCS Y1	*U803-49	Multi-site Y型多功能/1.1m	C4			
		LNOP DC-12	*U410-15	Adult finger 成人指夹/3.0m	C8			
PassportII, Exper,t Pssport XG, Accutor plus (*Non-Oximax)		LNOP DC-12	*U410S-15	Adult soft 成人软指套/3.0m	D1	No need 无需转接线		
			*U810-15	Multi-site Y型多功能/3.0m	C8			
		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
Accutor-Trio (*Non-Oximax)		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3	U708-69		
Trio Nouvelle Generation (* Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710-70		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
Digicare								
Digimax5500 (*BCI)		3044	U403-06	Adult finger 成人指夹/1.1m	C3	U708-76		
			U403S-06	Adult soft 成人软指套/1.1m	C5			
		3043	U803-06	Multi-site Y型多功能/1.1m	C3			
			U410-76	Adult finger 成人指夹/3.0m	C7			
			U410S-76	Adult soft 成人软指套/3.0m	C8			
			U810-76	Multi-site Y型多功能/3.0m	C7			
			U410-52	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-52	Adult soft 成人软指套/3.0m	C8			
			U810-52	Multi-site Y型多功能/3.0m	C7			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
Dima (Italia)								
(*Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
Dixtal(Brazil)-Philips								
Full compatible with Novamatrix Sensas		DX-1NC01-0	U410-113	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-113	Adult soft 成人软指套/3.0m	C8			
			U810-113	Multi-site Y型多功能/3.0m	C7			
Dolphin								
(*Digital Tech)			U403-50D	Adult finger 成人指夹/1.1m	C8	U708-50		
			U403S-50D	Adult soft 成人软指套/1.1m	D1			
			U803-50D	Multi-site Y型多功能/1.1m	C8			
(*Digital Tech)			U410-50D	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			U410S-50D	Adult soft 成人软指套/3.0m	D2			
			U810-50D	Multi-site Y型多功能/3.0m	D1			
Drager/Siemens德尔格/西门子								
Oxipac 2000 (*BCI)		3044	U403-06	Adult finger 成人指夹/1.1m	C3	No need 无需转接线		
			U403S-06	Adult soft 成人软指套/1.1m	C5			
		3043	U803-06	Multi-site Y型多功能/1.1m	C3			
Vitalert 1000,icero (*Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-28		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
Narkomed-Oxysat 2 (*Non-Oximax)			See the above	See the above	See the above	U710-85		
			See the above	See the above	See the above			
Sirecust, Sirem Spo2 (*Non-Oximax)			See the above	See the above	See the above	U708-33		
			See the above	See the above	See the above			
SC6000, 7000,9000 (*Non-Oximax)			See the above	See the above	See the above	U704-23 U708-23		
			See the above	See the above	See the above			
Oxytrend (*M tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-15		
		1863/ LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6			
		2258/ LNCS Y1	*U803-49	Multi-site Y型多功能/1.1m	C4			
		LNOP DC-12	*U410-15	Adult finger 成人指夹/3.0m	C8			
		LNOP DC-12	*U410S-15	Adult soft 成人软指套/3.0m	D1			
		*U810-15	Multi-site Y型多功能/3.0m	C8				

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
Edan理邦 (China)								
M9,M99A,M8,M8A,M8B, M3,M3B (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U708X-191		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
EDAN Oximetry M3,M3A,M3B, M8,M8A,M8B M9,M9A,M50,M80 iM8,iM8A,iM8B iM9,iM9A,iM50,iM80 (*Non-Oximax)		MS3-109069	U410-29	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-29	Adult soft 成人软指套/3.0m	C8			
			U810-29	Multi-site Y型多功能/3.0m	C7			
			U403-01	Adult finger 成人指夹/1.1m	C3			
M3,M3B,M8,M8A,M8B M9,M9A,M50,M80 iM9,iM9A,iM50,iM80 (*Oximax)		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5	U708-29		
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
			U410-191P	Adult finger 成人指夹/3.0m	C8			
(*Digital tech)			U410S-191P	Adult soft 成人软指套/3.0m	D1	No need 无需转接线		
			U810-191P	Multi-site Y型多功能/3.0m	C8			
			U410-46D	Adult finger 成人指夹/3.0m	C8			
H100B (*Non Oximax)		MS-109079	U403-01	Adult finger 成人指夹/1.1m	C3	No need 无需转接线		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
H100N (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	No need 无需转接线		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
Fukuda Denshi福田								
DS7100-DS7300 (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710-70		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
GE Corometrics								
188,129,129F, 128,128F (*Non Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
120,510,511, 118,556 (*Non Oximax)		See the above	See the above	See the above	See the above	U708-25 (4033CAX)		
128, 129 (*Non Oximax)		See the above	See the above	See the above	See the above	U708-21		
			U410-21	Adult finger 成人指夹/3.0m	C8			
			U410S-21	Adult soft 成人软指套/3.0m	D1			
			U810-21	Multi-site Y型多功能/3.0m	C8			
SERIE 250 (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710X-21		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
		See Oximax Marquette	See Oximax Marquette	See Oximax Marquette	See GE			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
GE Critikon								
Dianmap 9700/9701 9702 PLUS, MPS ORTABLE,/COMPACT, (* Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-01 (EC-8)		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
Dianmap,ProSeries V1 Pro 1000 V2 (Nellcor Optional)		See the above	See the above	See the above	See the above	U708-69 (SCP-10)		
Dianmap,ProCare, ProSeries V2Pro1000 V3 (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710-70 (DOC-10)		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
			U403-01	Adult finger 成人指夹/1.1m	C3			
Dianmap 8700/8710/8720Pulse, Oxyshuttle11 et 2		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	Not Available 不能提供		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
			U410-10	Adult finger 成人指夹/3.0m	C8			
			U410S-10	Adult soft 成人软指套/3.0m	D1			
			U810-10	Multi-site Y型多功能/3.0m	C8			
Dianmap ProCare Series V2 Pro 1000 V3 (*M Tech)		1863/LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-15 (2027263-001)		
		1863/LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6			
		2258/LNCS Y1	*U803-49	Multi-site Y型多功能/1.1m	C4			
		LNOP DC-12	*U410-15	Adult finger 成人指夹/3.0m	C8			
		LNOP DC-12	*U410S-15	Adult soft 成人软指套/3.0m	D1			
			*U810-15	Multi-site Y型多功能/3.0m	C8			
GE Datex-Ohmeda								
AS/3,CS/3, AS/5 (Built in M module) CH/S,CH/RSCG/2S, CARDIOCAPi-II Capnomac Ultima Light monitor Satellite, Oscar (Datex tech*)		OXY-F-DB	U403-09	Adult finger 成人指夹/1.1m	C3	U708-09		
			U103-09	Pediatric finger 儿童指夹/1.1m	C5			
			U403S-09	Adult soft 成人软指套/1.1m	C5			
			U103S-09	Pediatric soft 儿童软指套/1.1m	C5			
			U203S-09	Infant soft 婴儿软指套/1.1m	C5			
			U303-09	Neonate wrap 新生儿捆绑/1.1m	C5			
		OXY-W-DB	U603-09	Adult/Infant wrap成人/新生儿捆绑/1.1m	C3			
			U603B-09	Pedia/Infant wrap儿童/婴儿捆绑/1.1m	C3			
			U803-09	Multi-site Y型多功能/1.1m	C3			
		OXY-E-DB	U903-09	Adult ear clip 成人耳夹/1.1m	C5			

GE Datex-Ohmeda 下一页续接

GE Datex-Ohmeda 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头					
GE Datex-Ohmeda													
AS/3,CS/3, AS/5 (Built in M module) CH/S,CH/RSCG/2S, CARDIOQAPI-II Capnomac Ultima Light monitor Satellite, Oscar (Datex tech*)		OXY-F4-N	U410-09	Adult finger 成人指夹/3.0m	C8	No need 无需转接线							
			U110-09	Pediatric finger 儿童指夹/3.0m	D1								
			U410S-09	Adult soft 成人软指套/3.0m	D1								
			U110S-09	Pediatric soft 儿童软指套/3.0m	D1								
			U210S-09	Infant soft 婴儿软指套/3.0m	D1								
			U310-09	Neonate wrap 新生儿捆绑/3.0m	D1								
			OXY-W4-N	U610-09	Adult/Neo wrap 成人/新生儿捆绑/3.0m				C8				
				U610B-09	Pedia/infant wrap 儿童/婴儿捆绑/3.0m				C8				
			OXY-W4-N	U810-09	Multi-site Y型多功能/3.0m				C8				
			OXY-E4-N	U910-09	Adult ear clip 成人耳夹/3.0m				D1				
M-NSAT Module (AS/3,CS/3,S/5) (* Non Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-69 (MC-10)							
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5								
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3								
M-NSAT Module (4) (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710-70 (DOC-10)							
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7								
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5								
3700,3710,3770, 3775,Biox3740,3760, 3800,3900,3900P, 4700,5250RGB, Cardiocap/5, S/5, AS/3 (built in Module N-XOSAT) CS/3 (Built in module M-OSAT) (*Ohmeda)		OXY-F4-H	U410-02	Adult finger 成人指夹/3.0m	C8	No need 无需转接线							
			U110-02	Pediatric finger 儿童指夹/3.0m	D1								
			U410S-02	Adult soft 成人软指套/3.0m	D1								
			U110S-02	Pediatric soft 儿童软指套/3.0m	D1								
			U210S-02	Infant soft 婴儿软指套/3.0m	D1								
			U310-02	Neonate wrap 新生儿捆绑/3.0m	D1								
			OXY-W4-H	U610-02	Adult/Neo wrap 成人/新生儿捆绑/3.0m				C8				
			OXY-W4-H	U610B-02	Pedia/infant wrap 儿童/婴儿捆绑/3.0m				C8				
			OXY-W4-H	U810-02	Multi-site Y型多功能/3.0m				C8				
			OXY-E4-H	U910-02	Adult ear clip 成人耳夹/3.0m				D1				
		Ohmeda Tuffsat (*Ohmeda)		OXY-F1-H	U403-02				Adult finger 成人指夹/1.1m	C8	No need 无需转接线		
				OXY-F1-H	U403S-02				Adult soft 成人软指套/1.1m	D1			
		Ohmeda TruSat Ohmeda TruSat TD (*Oxyltip+ Technology)		OXY-F4-MC	U410-60				Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
	U110-60			Pediatric finger 儿童指夹/3.0m	D2								
	U410S-60			Adult soft 成人软指套/3.0m	D2								
	U110S-60			Pediatric soft 儿童软指套/3.0m	D2								
	U210S-60			Infant soft 婴儿软指套/3.0m	D2								
	U310-60			Neonate wrap 新生儿捆绑/3.0m	D2								
	OXY-W4-MC			U610-60	Adult/Neo wrap 成人/新生儿捆绑/3.0m	D1							
				U610B-60	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	D1							
				U810-60	Multi-site Y型多功能/3.0m	D1							
	OXY-E4-MC			U910-60	Adult ear clip 成人耳夹/3.0m	D2							

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
GE Hellige海力格								
Cardioserv Eagle 1000-3000 (*Non Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-55 (30344358)		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
Cardioserv (* Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	NA (2026751-001)		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
SMU,SMK,SMV (*Ohmeda)		See above *Ohmeda	See above *Ohmeda	See above *Ohmeda	See above *Ohmeda	No need 无需转接线		
GE Marquette								
Eagle 1000/3000/4000 Dash 2000/3000/4000 Solar 8000/8000M i9500(Tram 451/451N /851/851N) CardioLab/Mac-Lab, Corometrics 120/2120 (Tram 451/451N) Responder3000 (nelloor optional) (*Non Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-21 (2006644-001)		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
			U410-21	Adult finger 成人指夹/3.0m	C8			
			U110-21	Pediatric finger 儿童指夹/3.0m	D1			
			U410S-21	Adult soft 成人软指套/3.0m	D1			
			U110S-21	Pediatric soft 儿童软指套/3.0m	D1			
			U210S-21	Infant soft 婴儿软指套/3.0m	D1			
			U310-21	Neonate wrap 新生儿捆绑/3.0m	D1			
			U610-21	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8			
	U610B-21	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8					
	U810-21	Multi-site Y型多功能/3.0m	C8					
	U910-21	Adult ear clip 成人耳夹/3.0m	D1					
Dash2500/Dash 3000/ 4000/5000 Solar 8000/8000M i9500 (Tram 451N5/851N5) Corometrics 250 (GE-Oximax*) Ending " X" can be replaced by "P" for fully compatible with all Oximax oximetry		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710X-21 (2021406-001)		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
			U410-21P	Adult finger 成人指夹/3.0m	D1			
			U110-21P	Pediatric finger 儿童指夹/3.0m	D2			
			U410S-21P	Adult soft 成人软指套/3.0m	D2			
			U110S-21P	Pediatric soft 儿童软指套/3.0m	D2			
			U210S-21P	Infant soft 婴儿软指套/3.0m	D2			
			U310-21P	Neonate wrap 新生儿捆绑/3.0m	D2			
			U610-21P	Adult/Neo wrap 成人/新生儿捆绑/3.0m	D1			
	U610B-21P	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	D1					
	U810-21P	Multi-site Y型多功能/3.0m	D1					
	U910-21P	Adult ear clip 成人耳夹/3.0m	D2					
7000/7010,Tram 200/400/600/800 (*Ohmeda Tech.)		See above *Ohmeda	See above *Ohmeda	See above *Ohmeda	See above *Ohmeda			

GE Marquette 下页续接

GE Marquette 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
GE Marquette									
Dash 3000, 4000,5000 Solar 8000/8000M i9500 (Tram 451M,851M) Cardiolab/Mac-Lab (Tram 451M) Corometrics 259 (*M tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-21 (2027263-002)			
		1863/ LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6				
		2258/ LNCS Y1	*U803-49	Multi-site Y型多功能/1.1m	C4				
		*U410-21M	Adult finger 成人指夹/3.0m	C8	No need 无需转接线				
		*U110-21M	Pediatric finger 儿童指夹/3.0m	D1					
		*U410S-21M	Adult soft 成人软指套/3.0m	D1					
		*U110S-21M	Pediatric soft 儿童软指套/3.0m	D1					
		*U210S-21M	Infant soft 婴儿软指套/3.0m	D1					
		*U310-21M	Neonate wrap 新生儿捆绑/3.0m	D1					
		*U610-21M	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8					
		*U610B-21M	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8					
*U810-21M	Multi-site Y型多功能/3.0m	C8							
*U910-21M	Adult ear clip 成人耳夹/3.0m	D1							
GE Ohmeda OxyTIP+									
GE S/5 (4 E series module built in)		OXY-F4-GE	U410-121	Adult finger 成人指夹/3.0m	D1	No need 无需转接线			
		U110-121	Pediatric finger 儿童指夹/3.0m	D2					
		U410S-121	Adult soft 成人软指套/3.0m	D2					
		U110S-121	Pediatric soft 儿童软指套/3.0m	D2					
		U210S-121	Infant soft 婴儿软指套/3.0m	D2					
		U310-121	Neonate wrap 新生儿捆绑/3.0m	D2					
		U610-121	Adult/Neo wrap 成人/新生儿捆绑/3.0m	D1					
		U610B-121	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	D1					
		U810-121	Multi-site Y型多功能/3.0m	D1					
		U910-121	Adult ear clip 成人耳夹/3.0m	D2					
Ohmeda TruSat 3755,37XX series 3800,3900,3900P 4700, 5250 Cardiocap/5, S/5,AS/3 (N-XOSAT or M built in) CS/3 (M-OSAT built in) S/5 light, Satellite Puls II Capnomac UltimaS/5 , Coro 259 (*TruSignal tech)		OXY-F-UN	U403-127	Adult finger 成人指夹/1.1m	C8	Cables required depending on the monitor models 取决于机器端接头 参考页面血氧转接线			
		U103-127	Pediatric finger 儿童指夹/1.1m	D1					
		OXY-F-UN	U403S-127	Adult soft 成人软指套/1.1m	D1				
		U103S-127	Pediatric soft 儿童软指套/1.1m	D1					
		U203S-127	Infant soft 婴儿软指套/1.1m	D1					
		U303-127	Neonate wrap 新生儿捆绑/1.1m	D1					
		OXY-W-UN	U603-127	Adult/Neo wrap 成人/新生儿捆绑/1.1m	C8				
		OXY-W-UN	U603B-127	Pedia/infant wrap 儿童/婴儿捆绑/1.1m	C8				
		OXY-W-UN	U803-127	Multi-site Y型多功能/1.1m	C8				
		OXY-E-UN	U903-127	Adult ear clip 成人耳夹/1.1m	D1				

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
GE Trusignal								
		TS-F-D	U403-117	Adult finger 成人指夹/1.1m	C8	Required 必须配用转接线		
		U103-117	Pediatric finger 儿童指夹/1.1m	D1				
		TS-SA-D	U403S-117	Adult soft 成人软指套/1.1m	D1			
		U103S-117	Pediatric soft 儿童软指套/1.1m	D1				
		U203S-117	Infant soft 婴儿软指套/1.1m	D1				
		U303-117	Neonate wrap 新生儿捆绑/1.1m	D1				
		TS-W-D	U603-117	Adult/Neo wrap 成人/新生儿捆绑/1.1m	C8			
		TS-W-D	U603B-117	Pedia/infant wrap 儿童/婴儿捆绑/1.1m	C8			
		TS-W-D	U803-117	Multi-site Y型多功能/1.1m	C8			
		TS-E-D	U903-117	Adult ear clip 成人耳夹/1.1m	D1			
General 杰纳瑞(China)								
G1B (*Digital Tech)			U403-42D	Adult finger 成人指夹/1.1m	C5	No need 无需转接线		
			U403S-42D	Adult soft 成人软指套/1.1m	C7			
		G2A, G3A, G3G, G3D, G3F, G3G, G3H, G3L, G6A, G6A plus, G9L (*Digital Tech)	U410-65D	Adult finger 成人指夹/3.0m	C8			
		U410S-65D	Adult soft 成人软指套/3.0m	D1				
		U810-65D	Multi-site Y型多功能/3.0m	C8				
Goldway 金科威(China)								
UT4000A, UT4000B, UT4000F UT4000ACS			U410-37	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-37	Adult soft 成人软指套/3.0m	C8			
			U810-37	Multi-site Y型多功能/3.0m	C7			
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-37		
		See Pace tech	See Pace tech	See Pace tech	See Pace tech	No need 无需转接线		
UT4000A, UT4000C, UT4000E, UT4000F, UT4000APro, UT4000FPro, UT6000A			U410-11	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-11	Adult soft 成人软指套/3.0m	D1			
			U810-11	Multi-site Y型多功能/3.0m	C8			
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-11		
		See Pace tech	See Pace tech	See Pace tech	See Pace tech	No need 无需转接线		
3100		Digital sensor 数字探头	U410-61D	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
		U410S-61D	Adult soft 成人软指套/3.0m	D1				
		U810-67D	Multi-site Y型多功能/3.0m	C8				
UT4000B, UT4000F			U410-67D	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-67D	Adult soft 成人软指套/3.0m	D1			
			U810-67D	Multi-site Y型多功能/3.0m	C8			

Goldway 金科威(China) 下一页续接

Goldway金科威(China) 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
Goldway金科威(China)									
(*Oximax)			U410-38P	Adult finger 成人指夹/3.0m	D1	No need 无需转接线			
			U410S-38P	Adult soft 成人软指套/3.0m	D2				
			U810-38P	Multi-site Y型多功能/3.0m	D1				
UT6000A G3 (*Philips tech)			U410-91	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-91	Adult soft 成人软指套/3.0m	D1				
			U810-91	Multi-site Y型多功能/3.0m	C8				
Guoteng 国腾(China)									
PM300		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	No need 无需转接线			
			U403S-01	Adult soft 成人软指套/1.1m	C5				
			U410-171	Adult finger 成人指夹/3.0m	C8				No need 无需转接线
			U410S-171	Adult soft 成人软指套/3.0m	D1				
GT3000,GT5000, GT6000GT7000, GT9000,GT9002			U410-179	Adult finger 成人指夹/3.0m	C7	No need 无需转接线			
			U410S-179	Adult soft 成人软指套/3.0m	C8				
			U810-179	Multi-site Y型多功能/3.0m	C7				
IVY									
405A,405D (* Non Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U710-85			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5				
405P,405T (*Non Oximax)		See the above	See the above	See the above	See the above	U708-01			
		D-Y5	U803-01	Multi-site Y型多功能/1.1m	C3				
Vital-Guard 450CN(M) (*M tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-15			
			*U403S-49	Adult soft 成人软指套/1.1m	C6				
			1863/ LNCS DCI	*U803-49	Multi-site Y型多功能/1.1m				C4
		2258/ LNCS Y1	*U410-15	Adult finger 成人指夹/3.0m	C8				
		LNOP DC-12	*U410S-15	Adult soft 成人软指套/3.0m	D1				
	LNOP DC-12	*U810-15	Multi-site Y型多功能/3.0m	C8					
Kontron康强									
7840,7845,Minimon 7337Micromon, Supermon,Kolormon, Trakmon (*Non Oximax)			U410-13	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-13	Adult soft 成人软指套/3.0m	D1				
			U810-13	Multi-site Y型多功能/3.0m	C8				
		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3				
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5	U708-13			
		D-Y5	U803-01	Multi-site Y型多功能/1.1m	C3				
7150,7850, Ivy 2000(M)		See the above	See the above	See the above	See the above	U708-01			
		See the above	See the above	See the above	See the above	U708-53			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头					
Lohmeier隆美尔													
			U410-64	Adult finger 成人指夹/3.0m	D1	No need 无需转接线							
			U410S-64	Adult soft 成人软指套/3.0m	D2								
			U810-64	Multi-site Y型多功能/3.0m	D1								
*M Technology													
Rad-8 Rad-5/Rad-5v SatShare		1269/ LNOP DCI	*U403-125	Adult finger 成人指夹/1.1m	C7	Required 必须配用转接线							
		1276/ LNOP DCIP	*U103-125	Pediatric finger 儿童指夹/1.1m	C8								
		1269/ LNOP DCI	*U403S-125	Adult soft 成人软指套/1.1m	C8								
		1276/ LNOP DCIP	*U103S-125	Pediatric soft 儿童软指套/1.1m	C8								
			*U203S-125	Infant soft 婴儿软指套/1.1m	C8								
			*U303-125	Neonate wrap 新生儿捆绑/1.1m	C8								
			*U603-125	Adult/Neo wrap成人/新生儿捆绑/1.1m	C7								
			*U603B-125	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C7								
		1504/ LNOP Y1	*U803-125	Multi-site Y型多功能/1.1m	C7								
		LNOP TC-I	*U903-125	Adult ear clip 成人耳夹/1.1m	C8								
		1863/ LNCS DCI 1864/ LNCS DCIP 1863/ LNCS DCI			*U403-49				Adult finger 成人指夹/1.1m	C4	Required 必须配用转接线		
					*U103-49				Pediatric finger 儿童指夹/1.1m	C6			
					*U403S-49				Adult soft 成人软指套/1.1m	C6			
					*U103S-49				Pediatric soft 儿童软指套/1.1m	C6			
	*U203S-49			Infant soft 婴儿软指套/1.1m	C6								
	*U303-49			Neonate wrap 新生儿捆绑/1.1m	C6								
	*U603-49			Adult/Neo wrap成人/新生儿捆绑/1.1m	C4								
	*U603B-49			Pedia/infant wrap儿童/婴儿捆绑/1.1m	C4								
2258/ LNCS Y1	*U803-49			Multi-site Y型多功能/1.1m	C4								
LNCS TC-I	*U903-49			Adult ear clip 成人耳夹/1.1m	C6								
M Tech *Refurbished Only		LNOP DCSC	*U402-15	Adult finger 成人指夹/0.9m	C8	No need 无需转接线		Handheld					
		LNOP DC-12	*U410-15	Adult finger 成人指夹/3.0m	C8								
			*U110-15	Pediatric finger 儿童指夹/3.0m	D1								
		LNOP DC-12	*U410S-15	Adult soft 成人软指套/3.0m	D1								
			*U110S-15	Pediatric soft 儿童软指套/3.0m	D1								
			*U210S-15	Infant soft 婴儿软指套/3.0m	D1								
			*U310-15	Neonate wrap 新生儿捆绑/3.0m	D1								
			*U610-15	Adult/Neo wrap成人/新生儿捆绑/3.0m	C8								
	*U610B-15	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8										
	*U810-15	Multi-site Y型多功能/3.0m	C8										
	*U910-15	Adult ear clip 成人耳夹/3.0m	D1										

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
MDE/Invivo								
Escort 1 Escort 300 (*Non Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U710-85		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
20300,20401,20403, Prism (*Non Oximax)		See the above	See the above	See the above	See the above	U708-01		
Millenia, Escort,(M) M10TM,M12TM, M8TM (*M Tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-15		
		1863/ LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6			
		2258/ LNCS Y1	*U803-49	Multi-site Y型多功能/1.1m	C4			
		LNOP DC-12	*U410-15	Adult finger 成人指夹/3.0m	C8			
4500, ,Omni Track Omeg 1445Millennia			U410-14	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-14	Adult soft 成人软指套/3.0m	D1			
			U810-14	Multi-site Y型多功能/3.0m	C8			
4500, ,Omni Track Omeg 1445Millennia M Tech			U410-14M	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-14M	Adult soft 成人软指套/3.0m	D1			
			U810-14M	Multi-site Y型多功能/3.0m	C8			
Meditronic - Physio Control								
Lifepack 12 (*non Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
Lifepack 12, 20(M) (*M Tech)			*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-15		
			*U403S-49	Adult soft 成人软指套/1.1m	C6			
			*U803-49	Multi-site Y型多功能/1.1m	C4	No need 无需转接线		
			*U410-15	Adult finger 成人指夹/3.0m	C8			
			*U410S-15	Adult soft 成人软指套/3.0m	D1			
			*U810-15	Multi-site Y型多功能/3.0m	C8			
Medlab								
NANOx 10 POX 10 Perarl 10 (*BCI)		3044	U403-06	Adult finger 成人指夹/1.1m	C3	U708-76		
			U403S-06	Adult soft 成人软指套/1.1m	C5			
		3043	U803-06	Multi-site Y型多功能/1.1m	C3			
MEK								
MP 100,110, 400,1000			U410-26	Adult finger 成人指夹/3.0m	C6	No need 无需转接线		
			U410S-26	Adult soft 成人软指套/3.0m	C7			
			U810-26	Multi-site Y型多功能/3.0m	C6			
(*non Oximax)			U403-01	Adult finger 成人指夹/1.1m	C3	U708-56		
			U403S-01	Adult soft 成人软指套/1.1m	C5			
			U803-01	Multi-site Y型多功能/1.1m	C3			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
MEK									
			U410-56	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-56	Adult soft 成人软指套/3.0m	D1				
			U810-56	Multi-site Y型多功能/3.0m	C8				
Mennen 美能									
Envoy (*non Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-54			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5				
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3				
Horizon 1100,2000XL Mercury,MR1300 (*non Oximax)		See the above	See the above	See the above	See the above	U710-85			
Horizon XL,Mercury 1200,MR1300 (*non Oximax)		See the above	See the above	See the above	See the above	U708-24			
Horizon XL (*Ohmeda)			OXY-F4-H	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-02	Adult soft 成人软指套/3.0m	D1				
			U810-02	Multi-site Y型多功能/3.0m	C8				
Horizon,Envoy, Vitalogik (M) (*M Tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-59			
			*U403S-49	Adult soft 成人软指套/1.1m	C6				
		2258/ LNCS Y1	*U803-49	Multi-site Y型多功能/1.1m	C6				
Mindray迈瑞 (China)									
509B, MEC-1000, MEC-1200 New PM6000 PM7000/8000/9000, VS 800,PM600 PM9000 Express (*Mindray Tech)			U410-29	Adult finger 成人指夹/3.0m	C7	No need 无需转接线			
			U110-29	Pediatric finger 儿童指夹/3.0m	C8				
			U410S-29	Adult soft 成人软指套/3.0m	C8				
			U110S-29	Pediatric soft 儿童软指套/3.0m	C8				
			U210S-29	Infant soft 婴儿软指套/3.0m	C8				
			U310-29	Neonate wrap 新生儿捆绑/3.0m	C8				
			U610-29	Adult/Neo wrap成人/新生儿捆绑/3.0m	C7				
			U610B-29	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C7				
			U810-29	Multi-site Y型多功能/3.0m	C7				
			U910-29	Adult ear clip成人耳夹/3.0m	C8				
			U403-01	Adult finger 成人指夹/1.1m	C3	U708-29 (0010-20-42594)			
			512D	U403S-01	Adult soft 成人软指套/1.1m				C5
			DS100A	U803-01	Multi-site Y型多功能/1.1m				C3
PM-5000,PM 6000 PM-8000 (*Mindray Tech)			See the above	See the above	See the above	U708-22 (0515-30-11221)			
			U410-22	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-22	Adult soft 成人软指套/3.0m	D1				
			U810-22	Multi-site Y型多功能/3.0m	C8				

Mindray迈瑞 (China) 下一页继续 >

← Mindray迈瑞 (China) 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头			
Mindray迈瑞 (China)											
PM-6000 (*M tech)		1863/LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-22 (0010-30-12452)					
		1863/LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6						
		2258/LNCS YI	*U803-49	Multi-site Y型多功能/1.1m	C4		No need 无需转接线				
			*U410-22M	Adult finger 成人指夹/3.0m	C8						
New PM-6000, PM-8000 (*M tech)			*U410S-22M	Adult soft 成人软指套/3.0m	D1	No need 无需转接线					
			*U810-22M	Multi-site Y型多功能/3.0m	C8						
			*U410-29M	Adult finger 成人指夹/3.0m	C7						
New PM-6000, PM7000,8000,9000 PM-9000 Express VS 800 New PM6000 (*M tech)		See the above M sensor	See the above M sensor	See the above M sensor	See the above	U708M-30					
			*U410-30M	Adult finger 成人指夹/3.0m	C7				No need 无需转接线		
			*U410S-30M	Adult soft 成人软指套/3.0m	C8						
New PM-6000 PM7000,8000,9000 VS 800 PM-9000 Express (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710X-180					
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7						
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5				No need 无需转接线		
			U410-180P	Adult finger 成人指夹/3.0m	C8						
Beneview T5/T8 (*Mindray Tech)		See the above Non Oximax Sensor	See the above Non Oximax Sensor	1.1m	See the above	U708-48					
			U410-48	Adult finger 成人指夹/3.0m	C8				No need 无需转接线		
			U410S-48	Adult soft 成人软指套/3.0m	D1						
			U810-48	Multi-site Y型多功能/3.0m	C8						
Beneview T5/T8 (Oximax*)		See Nellcor Oximax	See Nellcor Oximax	1.1m	See Oximax	U710X-48					
Beneview T5/T8 (*M tech)		See M	See M	1.1m	See M	U708M-48					
503,507,508, 509,PM9000 (*Criticare tech)			U410-05	Adult finger 成人指夹/3.0m	C8	No need 无需转接线					
			U410S-05	Adult soft 成人软指套/3.0m	D1						
			U810-05	Multi-site Y型多功能/3.0m	C8						
500A,509B, PM9000 (*BCI)			U410-31	Adult finger 成人指夹/3.0m	C8	No need 无需转接线					
			U410S-31	Adult soft 成人软指套/3.0m	D1						
			U810-31	Multi-site Y型多功能/3.0m	C8						
PM-50,PM-60 (*Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	No need 无需转接线					
	DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5							

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头			
Minolta美能达											
PULSOX-300 PULSOX-300i		SR-5C	U401-105	Adult finger 成人指夹/0.3m	C8	No need 无需转接线					
		SR-5C	U401S-105	Adult soft 成人软指套/0.3m	D1						
M&B麦邦(China)											
MB526,MB526+,MB526T8 MB526T12,MB526+hw MB800,CD2000,CD2000+			U410-92	Adult finger 成人指夹/3.0m	C7	No need 无需转接线					
			U410S-92	Adult soft 成人软指套/3.0m	C8						
			U810-92	Multi-site Y型多功能/3.0m	C7						
CD2000,CD2000+ (* Oximax tech)			U410-92P	Adult finger 成人指夹/3.0m	C8	No need 无需转接线					
			U410S-92P	Adult soft 成人软指套/3.0m	D1						
			U810-92P	Multi-site Y型多功能/3.0m	C8						
Nellcor-Tyco											
NPB 40, NPB20, NPB75 NPB 180, NPB 185, NPB-190, NPB-195, NPB 3900, NPB 4000 (*Non-Oximax)+		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U704-01 U708-01					
		OXI-P/I	U103-01	Pediatric finger 儿童指夹/1.1m	C5						
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5						
		OXI-P/I	U103S-01	Pediatric soft 儿童软指套/1.1m	C5						
		OXI-P/I	U203S-01	Infant soft 婴儿软指套/1.1m	C5						
		OXI-A/N	U303-01	Neonate wrap 新生儿捆绑/1.1m	C5						
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3						
		OXI-A/N	U603-01	Adult/Neo wrap成人/新生儿捆绑/1.1m	C3						
		OXI-P/I	U603B-01	Pedia/Infant wrap儿童/婴儿捆绑/1.1m	C3						
		D-YS+YSE	U903-01	Adult ear clip 成人耳夹/1.1m	C5						
			U410-01	Adult finger 成人指夹/3.0m	C5				No need 无需转接线		
			U110-01	Pediatric finger 儿童指夹/3.0m	C7						
	U410S-01	Adult soft 成人软指套/3.0m	C7								
	U110S-01	Pediatric soft 儿童软指套/3.0m	C7								
	U210S-01	Infant soft 婴儿软指套/3.0m	C7								
	U310-01	Neonate wrap 新生儿捆绑/3.0m	C7								
	U610-01	Adult/Neo wrap成人/新生儿捆绑/3.0m	C5								
	U610B-01	Pedia/Infant wrap儿童/婴儿捆绑/3.0m	C5								
	U810-01	Multi-site Y型多功能/3.0m	C5								
	U910-01	Adult ear clip 成人耳夹/3.0m	C7								
N-100,N-200		See the above	See the above	1.1m	See the above	U710-85					

Nellcor-Tyco 下一页续接 →














↳ Nellcor-Tyco续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
Nellcor-Tyco								
NPB3000-3100, Symphony, NPB-290, NPB-295, NPB-395		See the above	See the above	1.1m	See the above	U708-69		
			U410-69	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U110-69	Pediatric finger 儿童指夹/3.0m	D1			
			U410S-69	Adult soft 成人软指套/3.0m	D1			
			U110S-69	Pediatric soft 儿童软指套/3.0m	D1			
			U210S-69	Infant soft 婴儿软指套/3.0m	D1			
			U310-69	Neonate wrap 新生儿捆绑/3.0m	D1			
			U610-69	Adult/Neo wrap成人/新生儿捆绑/3.0m	C8			
			U610B-69	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8			
			U810-69	Multi-site Y型多功能/3.0m	C8			
	U910-69	Adult ear clip 成人耳夹/3.0m	D1					
N-65 (*Oximax) Ending "X" can be replaced by "XP" or "P" for fully compatible with all Oximax oximetry		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	No need 无需转接线		
		OXI-P/I	U103-01X	Pediatric finger 儿童指夹/1.1m	C7			
		OXI-A/N	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		OXI-P/I	U103S-01X	Pediatric soft 儿童软指套/1.1m	C7			
		OXI-P/I	U203S-01X	Infant soft 婴儿软指套/1.1m	C7			
		OXI-A/N	U303-01X	Neonate wrap 新生儿捆绑/1.1m	C7			
		OXI-A/N	U603-01X	Adult/Neo wrap成人/新生儿捆绑/1.1m	C5			
		OXI-P/I	U603B-01X	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C5			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
		D-YS+YSE	U903-01X	Adult ear clip 成人耳夹/1.1m	C7			
N550,N560,N595, N600,N751, N5500		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U701-01X		
		See Nellcor Oximax	See Nellcor Oximax	1.1m	See Oximax	U710-70		
Neusoft 东软(China)								
NSC series (*Oximax Tech)			U410-78P	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			U410S-78P	Adult soft 成人软指套/3.0m	D2			
			U810-78P	Multi-site Y型多功能/3.0m	D1			
			U410-107	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-107	Adult soft 成人软指套/3.0m	C8			
			U810-107	Multi-site Y型多功能/3.0m	C7			
			U810-107	Multi-site Y型多功能/3.0m	C7			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
New tech组泰克 (China)									
Pulse Oximer NT1A,NT1B,NT1D (Digital Sensor)			U403-42D	Adult finger 成人指夹/1.1m	C5	No need 无需转接线			
			U403S-42D	Adult soft 成人软指套/1.1m	C7				
Patient monitor Digital sensor			U410-65D	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-65D	Adult soft 成人软指套/3.0m	D1				
NT2A,NT2C			U810-65D	Multi-site Y型多功能/3.0m	C8	No need 无需转接线			
			U410-165D	Adult finger 成人指夹/3.0m	D1				
			U410S-165D	Adult soft 成人软指套/3.0m	D2				
			U810-165D	Multi-site Y型多功能/3.0m	D1				
NT3A,NT3B, NT3E,NT3F, NT3B60			U410-167	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-167	Adult soft 成人软指套/3.0m	D1				
			U810-167	Multi-site Y型多功能/3.0m	C8				
NT3A,NT3B, NT3E,NT3F,NT3B60 (*Digital Tech)			U410-173D	Adult finger 成人指夹/3.0m	D1	No need 无需转接线			
			U410S-173D	Adult soft 成人软指套/3.0m	D2				
			U810-173D	Multi-site Y型多功能/3.0m	D1				
Nihon Kohden光电									
BSM-1101/2101 /23101 7100,7200, OLV-1100 ALife scope			U410-16	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U110-16	Pediatric finger 儿童指夹/3.0m	D1				
			U410S-16	Adult soft 成人软指套/3.0m	D1				
			U110S-16	Pediatric soft 儿童软指套/3.0m	D1				
			U210S-16	Infant soft 婴儿软指套/3.0m	D1				
			U310-16	Neonate wrap 新生儿捆绑/3.0m	D1				
			U610-16	Adult/Neo wrap成人/新生儿捆绑/3.0m	D1				
			U610B-16	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8				
			U810-16	Multi-site Y型多功能/3.0m	C8				
			U910-16	Adult ear clip 成人耳夹/3.0m	D1				
Lifescop S, BSS-9800 Lifescop M, BSM 9510 Lifescop P BSM 4100 Life scope A, BSM 5100 Lifescop L, BSM 2350 Lifescop I BSM 2301K, 2303K		TL-201T	U405-16	Adult finger 成人指夹/1.5m	C7	U710-36			
			U105-16	Pediatric finger 儿童指夹/1.5m	C8				
			U405S-16	Adult soft 成人软指套/1.5m	C8				
			U105S-16	Pediatric soft 儿童软指套/1.5m	C8				
			U205S-16	Infant soft 婴儿软指套/1.5m	C8				
			U305-16	Neonate wrap 新生儿捆绑/1.5m	C8				
			U605-16	Adult/Neo wrap成人/新生儿捆绑/1.5m	C7				
			U605B-15	Pedia/infant wrap儿童/婴儿捆绑/1.5m	C7				
		TL-220T	U805-16	Multi-site Y型多功能/1.5m	C7				
			U905-16	Adult ear clip 成人耳夹/1.5m	C8				

Nihon Kohden光电 下一页续接 ↳

Nihon Kohden光电 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
Nihon Kohden光电								
Lifescope S, BSS-9800 Lifescope M, BSM 9510 Lifescope P BSM 4100 Life scope A, BSM 5100 Lifescope L, BSM 2350 Lifescope I BSM 2301K, 2303K			U410-36	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			U110-36	Pediatric finger 儿童指夹/3.0m	D2			
			U410S-36	Adult soft 成人软指套/3.0m	D2			
			U110S-36	Pediatric soft 儿童软指套/3.0m	D2			
			U210S-36	Infant soft 婴儿软指套/3.0m	D2			
			U310-36	Neonate wrap 新生儿捆绑/3.0m	D2			
			U610-36	Adult/Neo wrap成人/新生儿捆绑/3.0m	D1			
			U610B-36	Pedia/infant wrap儿童/婴儿捆绑/3.0m	D1			
			U810-36	Multi-site Y型多功能/3.0m	D1			
	U910-36	Adult ear clip 成人耳夹/3.0m	D2					
OPV-1500		The same as above	The same as above	1.1m	The same as above	U710-106		
Nonin燕牌								
8500 2500 9840 Avant 4000,9600, 9700 Avant 2120 PulseSense Model 7500,7500FO (Nonin*)		8000AA	U403-08	Adult finger 成人指夹/1.1m	C3	No need 无需转接线		
		8000AP	U103-08	Pediatric finger 儿童指夹/1.1m	C5			
		8000SL	U403S-08	Adult soft 成人软指套/1.1m	C5			
		8000SM	U103S-08	Pediatric soft 儿童软指套/1.1m	C5			
			U203S-08	Infant soft 婴儿软指套/1.1m	C5			
			U303-08	Neonate wrap 新生儿捆绑/1.1m	C5			
			U603-08	Adult/Neo wrap成人/新生儿捆绑/1.1m	C5			
			U603B-08	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C3			
			U803-08	Multi-site Y型多功能/1.1m	C3			
			U903-08	Adult ear clip 成人耳夹/1.1m	C5			
			U410-08L	Adult finger 成人指夹/3.0m	C5			
8600		The same as above	The same as above	1.1m	The same as above	NA		
8604		The same as above	The same as above	1.1m	The same as above	U708-08		
				U410-08	Adult finger 成人指夹/3.0m	C8	No need 无需转接线	
Avant 4000 Wearable Digital Pulse Oximetry			U401-146	Adult finger 成人指夹/0.3m	D1	No need 无需转接线		
			U101-146	Pediatric finger 儿童指夹/0.3m	D1			
			U401S-146	Adult soft 成人软指套/0.3m	D1			
			U101S-146	Pediatric soft 儿童软指套/0.3m	D1			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
Novamatrix诺医								
512/513/710/715, MARS			U403-73	Adult finger 成人指夹/1.1m	C5	No need 无需转接线		
			U103-73	Pediatric finger 儿童指夹/1.1m	C7			
			U403S-73	Adult soft 成人软指套/1.1m	C7			
			U103S-73	Pediatric soft 儿童软指套/1.1m	C7			
			U203S-73	Infant soft 婴儿软指套/1.1m	C7			
			U303-73	Neonate wrap 新生儿捆绑/1.1m	C7			
			U603-73	Adult/Neo wrap成人/新生儿捆绑/1.1m	C5			
			U603B-73	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C5			
			U803-73	Multi-site Y型多功能/1.1m	C5			
			U903-73	Adult ear clip 成人耳夹/1.1m	C7			
			U403-103	Adult finger 成人指夹/1.1m	C8			
	U103-103	Pediatric finger 儿童指夹/1.1m	D1					
	U403S-103	Adult soft 成人软指套/1.1m	D1					
	U103S-103	Pediatric soft 儿童软指套/1.1m	D1					
	U203S-103	Infant soft 婴儿软指套/1.1m	D1					
	U303-103	Neonate wrap 新生儿捆绑/1.1m	D1					
	U603-103	Adult/Neo wrap成人/新生儿捆绑/1.1m	C8					
	U603B-103	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C8					
	U803-103	Multi-site Y型多功能/1.1m	C8					
	U903-103	Adult ear clip 成人耳夹/1.1m	D1					
505/510/511 /515/B/C/520 7000/7100,Oxypleth (Novamatrix*)			U410-03	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U110-03	Pediatric finger 儿童指夹/3.0m	D1			
			U410S-03	Adult soft 成人软指套/3.0m	D1			
			U110S-03	Pediatric soft 儿童软指套/3.0m	D1			
			U210S-03	Infant soft 婴儿软指套/3.0m	D1			
			U310-03	Neonate wrap 新生儿捆绑/3.0m	D1			
			U610-03	Adult/Neo wrap成人/新生儿捆绑/3.0m	C8			
			U610B-03	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8			
			U810-03	Multi-site Y型多功能/3.0m	C8			
			U910-03	Adult ear clip 成人耳夹/3.0m	D1			
PaceTech/Generra								
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-38		
			U410-38	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
		U410S-38	Adult soft 成人软指套/3.0m	C8				
		U810-38	Multi-site Y型多功能/3.0m	C7				

PaceTech/Generra 下一页续接

↳ PaceTech/Generra 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
PaceTech/Generra								
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-39		
			U410-39	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-39	Adult soft 成人软指套/3.0m	D1			
			U810-39	Multi-site Y型多功能/3.0m	C8			
Mediaid/Palco								
M300/M305/340, M5305/M5340			U410-18	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U110-18	Pediatric finger 儿童指夹/3.0m	D1			
			U410S-18	Adult soft 成人软指套/3.0m	D1			
			U110S-18	Pediatric soft 儿童软指套/3.0m	D1			
			U210S-18	Infant soft 婴儿软指套/3.0m	D1			
			U310-18	Neonate wrap 新生儿捆绑/3.0m	D1			
			U610-18	Adult/Neo wrap成人/新生儿捆绑/3.0m	C8			
			U610B-18	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8			
			U810-18	Multi-site Y型多功能/3.0m	D1			
M30/40,M100,M130, M900,M960P, M6 Aster,M4 Nova, M1000			U410-148	Adult finger 成人指夹/3.0m	D1.3	No need 无需转接线		
			U110-148	Pediatric finger 儿童指夹/3.0m	D1.6			
			U410S-148	Adult soft 成人软指套/3.0m	D1.6			
			U110S-148	Pediatric soft 儿童软指套/3.0m	D1.6			
			U210S-148	Infant soft 婴儿软指套/3.0m	D1.6			
			U310-148	Neonate wrap 新生儿捆绑/3.0m	D1.6			
			U610-148	Adult/Neo wrap成人/新生儿捆绑/3.0m	D1.3			
			U610B-148	Pedia/infant wrap儿童/婴儿捆绑/3.0m	D1.3			
			U810-148	Multi-site Y型多功能/3.0m	D1.6			
Philips/HP飞利浦								
78352C,78354C, 78834C,M1020A, M1025B,M1025A, M1350B,M1772A/B, M1732A/B,M2475B (*Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-20		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-Y5	U803-01	Multi-site Y型多功能/1.1m	C3	U708-40		
		M1196A	U403-91 U410-91	Adult finger 成人指夹/1.1m 3.0m	C8			
		M1192A	U103-91	Pediatric finger 儿童指夹/1.1m	D1			
		M1191A	U403S-91 U410S-91	Adult soft 成人软指套/1.1m 3.0m	D1			
		M1192A	U103S-91	Pediatric soft 儿童软指套/1.1m	D1			
		M1195A	U203S-91	Infant soft 婴儿软指套/1.1m	D1			
		M1193A	U303-91	Neonate wrap 新生儿捆绑/1.1m	D1			
			U603-91	Adult/Neo wrap成人/新生儿捆绑/1.1m	C8			
			U603B-91	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C8			
			U803-91	Multi-site Y型多功能/1.1m	D1			
		M1194A	U903-91	Adult ear clip 成人耳夹/1.1m	D1			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
Philips/HP飞利浦								
78352C,78354C, 78834C,M1020A, M1025B,M1025A, M1350B,M1772A/B, M1732A/B,M2475B (*Non-Oximax)		M1190A	U410-20	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U110-20	Pediatric finger 儿童指夹/3.0m	D1			
		M1190A	U410S-20	Adult soft 成人软指套/3.0m	D1			
			U110S-20	Pediatric soft 儿童软指套/3.0m	D1			
			U210S-20	Infant soft 婴儿软指套/3.0m	D1			
			U310-20	Neonate wrap 新生儿捆绑/3.0m	D1			
			U610-20	Adult/Neo wrap成人/新生儿捆绑/3.0m	C8			
			U610B-20	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8			
			U810-20	Multi-site Y型多功能/3.0m	C8			
			U910-20	Adult ear clip 成人耳夹/3.0m	D1			
C1,C3,M3,M4, VirDia, MP20,MP30, MP40/50/60/70/90, VM4,VM6,VM8 (Philips*)		See Philips' sensor	See Philips' sensor	See Philips' sensor	See Philips' sensor	U708-41		
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-43		
A1, A3 (*Non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m/3.0m	See Nellcor	U708-01		
MP20/30/40/50 PM60/70/90 VN4,VM6,VM8 (*Oximax)		DS100A	U403-01P	Adult finger 成人指夹/1.1m	C5	U710X-43		
		DS100A	U403S-01P	Adult soft 成人软指套/1.1m	C7			
		D-Y5	U803-01P	Multi-site Y型多功能/1.1m	C5			
MP20/30/40/ PM50/60/70/90 VN4,VM6,VM8 (*M tech)		See M	See M	1.1m	See M	U708M-43		
			*U410-91M	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			*U410S-91M	Adult soft 成人软指套/3.0m	D2			
			*U810-91M	Multi-site Y型多功能/3.0m	D1			
RGB medical								
(*BCI)		3044	U403-06	Adult finger 成人指夹/1.1m	C3	U708-76		
			U403S-06	Adult soft 成人软指套/1.1m	C3			
		3078	U803-06	Adult ear clip 成人耳夹/1.1m	C5			
Schiller ODM/Bruker								
CM8 SM783,784,785, TM910		See RGB medical	See RGB medical	See RGB medical	See RGB medical	U708-17		
			U410-17	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-17	Adult soft 成人软指套/3.0m	D1			
			U810-17	Multi-site Y型多功能/3.0m	C8			

Schiller ODM/Bruker 下页续接

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
Schiller ODM/Bruker									
TM7, Argus LCM, (*M tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-17			
		1863/ LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6				
		2258/ LNCS Y1	*U803-49	Multi-site Y型多功能/1.1m	C4				
			*U410-17M	Adult finger 成人指夹/3.0m	C8				
			*U410S-17M	Adult soft 成人软指套/3.0m	D1				
Propaq CS, Encore, Micropaq (New) (M) (*M tech)		See M	See M	1.1m	See M	U708M-15			
		LNOP DC-12	*U410-15	Adult finger 成人指夹/3.0m	C8				No need 无需转接线
		LNOP DC-12	*U410S-15	Adult soft 成人软指套/3.0m	D1				
			*U810-15	Multi-site Y型多功能/3.0m	C8				
Propaq CS, Encore Series, Atlas (* non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-01			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5				
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3				
TM 910-Defigard 500(M) (* non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m/3.0m	See Nellcor	U708-63			
Argus VCM (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710-70			
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7				
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5				
Spacelabs太空									
90419/01, 90651A/07, 90501, Ultraview 90496 (*Novamatrix)			U410-03	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U110-03	Pediatric finger 儿童指夹/3.0m	D1				
			U410S-03	Adult soft 成人软指套/3.0m	D1				
			U110S-03	Pediatric soft 儿童软指套/3.0m	D1				
			U210S-03	Infant soft 婴儿软指套/3.0m	D1				
			U310-03	Neonate wrap 新生儿捆绑/3.0m	D1				
			U610-03	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8				
			U610B-03	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8				
			U810-03	Multi-site Y型多功能/3.0m	C8				
			U910-03	Adult ear clip 成人耳夹/3.0m	D1				
90351-0/6, 90465, 90466, 90456, PC exoress (*Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-04			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5				
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3				
			U410-04	Adult finger 成人指夹/3.0m	C8				
			U110-04	Pediatric finger 儿童指夹/3.0m	D1				
			U410S-04	Adult soft 成人软指套/3.0m	D1				
			U110S-04	Pediatric soft 儿童软指套/3.0m	D1				
			U210S-04	Infant soft 婴儿软指套/3.0m	D1				
			U310-04	Neonate wrap 新生儿捆绑/3.0m	D1				
			U610-04	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8				
	U610B-04	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8						
	U810-04	Multi-site Y型多功能/3.0m	C8						
	U910-04	Adult ear clip 成人耳夹/3.0m	D1						

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
Spacelabs太空									
Ultraview (*Non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-74			
			U410-74	Adult finger 成人指夹/3.0m	D1				No need 无需转接线
			U410S-74	Adult soft 成人软指套/3.0m	D2				
			U810-74	Multi-site Y型多功能/3.0m	D1				
			U410-74M	Adult finger 成人指夹/3.0m	D1				
Ultraview (M Tech)			U410S-74M	Adult soft 成人软指套/3.0m	D2	No need 无需转接线			
			U810-74M	Multi-site Y型多功能/3.0m	D1				
			U410-74P	Adult finger 成人指夹/3.0m	D1				
Ultraview (Oximax Tech)			U410S-74P	Adult soft 成人软指套/3.0m	D2	No need 无需转接线			
			U810-74P	Multi-site Y型多功能/3.0m	D1				
			U403-06	Adult finger 成人指夹/1.1m	C3				
PRO2 (*BCI)		3044	U403S-06	Adult soft 成人软指套/1.1m	C5	No need 无需转接线			
		3043	U803-06	Multi-site Y型多功能/1.1m	C3				
Zoll									
M series AED (*M tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-15			
		1863/ LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6				
		2258/ LNCS Y1	*U803-49	Multi-site Y型多功能/1.1m	C4				
		LNOP DC-12	*U410-15	Adult finger 成人指夹/3.0m	C8				
		LNOP DC-12	*U410S-15	Adult soft 成人软指套/3.0m	D1				
		*U810-15	Multi-site Y型多功能/3.0m	C8					
Zondan和心重典 (China)									
Apollo N2, Apollo N4 Apollo N5 (M tech optional) DS120C			U410-45	Adult finger 成人指夹/3.0m	C7	No need 无需转接线			
			U410S-45	Adult soft 成人软指套/3.0m	C8				
			U810-45	Multi-site Y型多功能/3.0m	C7				
ZD120B ZD120D ZD120E (*Digital Tech)			U410-66D	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-66D	Adult soft 成人软指套/3.0m	D1				
			U810-66D	Multi-site Y型多功能/3.0m	C8				
Unicare(China)									
(*Digital tech)			U410-66D	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-66D	Adult soft 成人软指套/3.0m	D1				
			U810-66D	Multi-site Y型多功能/3.0m	C8				

Please be Reminded : The other types of sensors like U100, U100S, U200S, U300, U600, U600B and U900 of the listed respected brands are available upon demands.

注意：其他类型的兼容于各厂家的探头，如U100儿童软指套，U200婴儿软指套，U300新生儿捆绑式，U600包裹式及U900成人耳夹探头都可以选择。

