



INNOVATIVE
Medical Solutions
Specialists in Medical Electronics

SpO₂ Interface Cables

1 July 2015

Phone: 1300 884 612
Fax: 1300 884 613
Email: sales@innovativemedical.com.au



③ Interface Cable

血氧转接电缆

Specification 产品规格

Cable material 电缆材料	TPU Jacket, Light gray, Blue TPU 外壳, 浅灰色, 蓝色
Connector 插头	Nylon plug & Golden plated pins 尼龙 插头 & 镀金插针
Cable length 电缆长度	1.1m/2.2m/3.0m
Cable diameter 电缆直径	4.0mm/5.0mm
Package 包装	1 pcs per bag 每袋一包



REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
Artema 阿替玛						
U708-80	59831-0	Nellcor,non-Oximax	2.2m			C3
BCI						
U708-76		BCI	2.2m			B8
U708-06		BCI	2.2m			B8
Biolight 宝莱特						
U710X-51		Nellcor Oximax	3.0m			C3
U708-68		Digital sensor 数字探头	2.2m			C3
Bionet						
U708-58		Nellcor,non-Oximax	2.2m			B8
U708-108		Bionet	2.2m			C3
Criticare						
U708-05	518LD	CSI	2.2m			C3
U708-75	518DD	CSI	2.2m			B8
Datascope						
U708-07	0012-00-0516-02	DATASCOPE	2.2m			C3
Dolphin 海豚						
U708-50		Dolphin digital sensor Dolphin数字探头	2.2m			C3
Drager 德尔格						
U708-28	5720084	Nellcor,non-Oximax	2.2m			D4
EDAN 理邦						
U708X-191	MS1-30131 (M9,M9A,M8,M8A,M8B)	Nellcor-Oximax	2.2m			C3
U708-29	M13-30691 (M9,M9A,M8,M8A,M8B, M3,M3B)	EDAN Oximetry	2.2m			B8

REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
GE healthcare						
U708-21	2006644-001 E9004GE	Nellcor,non-Oximax	2.2M			C3
U710X-21	2021406-001	Nellcor,Oximax	3.0M			C7
*U708-21M	M technology	Nellcor Non-Oximax	2.2M			C3
*U708M-21	LNC-10-GE 2027263-002 (REF 2016)	M tech.(Masimo LNCS)	2.2M			C3
*U708M-21P	20169041-001	M tech.(Masimo LNOP)	2.2M			C7
U708-09	OXY-C3	Datex	2.2m			C3
U708-81	Datex>>Oxytip Used with Oxy-C3 cable	Oxytip	2.2m			C3
U708-55	30344358	Hellige(Nellcor)	2.2m			C3
U708-25	4033CAX	Corometrics(Nellcor)	2.2m			C3
U708-72*	All Ohmeda Oximeters	Ohmeda Used with Unimed N500 or U500 Disposable sensors	2.2m			C3
U708-160	OXY-MC3	GE OxyTip+ sensor	2.2m			C7
U708-102 U704-102	OXY-OL3	GE OxyTip+ sensor	2.2m 1.1m			C5
U708-209	OXY-SL3	GE OxyTip+ sensor	2.2m			C5
U708-121 U710-121	OXY-ES3	GE OxyTip+ sensor	2.2m 3.0m			C7
U708-81 U701-81	OXY-SLA OXY-SLC	GE OxyTip+ sensor	2.2m 0.3m			C5
U708-111	TS-G3	GE OxyTip+ sensor	2.2m			C5
U708-112	TS-H3	GE OxyTip+ sensor	2.2m			C5
U708-116	TS-M3	GE OxyTip+ sensor	2.2m			C7
U708-119	TS-N3	GE OxyTip+ sensor	2.2m			C5
Goldway 金科威						
U708-11		Nellcor Non-Oximax	2.2m			C3
U708-37		Nellcor Non-Oximax	2.2m			B8

REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
Goldway 金科威						
U708-38		Nellcor Non-Oximax	2.2m			B8
U710X-38		Nellcor Oximax	3.0m			C3
Invivo 鹰威						
U708-14		Nellcor,non-Oximax	2.2m			C3
*U708-14M		(Masimo) Nellcor, non-Oximax	2.2m			C3
*U708M-14		Masimo LNCS sensor	2.2m			C3
Kontron 康强						
U708-13		Nellcor,non-Oximax	2.2m			C3
U708-53		Nellcor,non-Oximax	2.2m			B8
M Technology						
*U708-15	M technology	(Masimo) Nellcor, non-Oximax	2.2m			C3
*U708M-15	LNCS Patient Cable LNC-4(REF 2017) LNC-10(REF 1814) LNC-14(REF 2013)	Masimo LNCS sensor	2.2m			C3
*U708M-15P	PC 04 / PC 08 / PC 12	Masimo LNOP sensor	2.2m			C3
*U704M-49	LNC-4-Ext (REF 2021) For LNCS Cable	Masimo LNCS sensor	1.1m			B8
*U701M-125	LNCS to LNOP LNCS to PC series (REF 1816)	Masimo LNCS sensor	0.3m			C3
MEK						
U708-56		MEK (Nellcor non-oximax)	2.2m			C3
Mennen 美能						
U708-54		Nellcor,non-Oximax	2.2m			C3
U708-24		Nellcor,non-Oximax	2.2m			D4
*U708-59M	M technology	M tech. (Nellcor, non-Oximax)	2.2m			C3
*U708M-59	M technology	Masimo LNCS sensor	2.2m			C3
U710X-59		Nellcor Oximax	3.0m			C3

REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
Mindray迈瑞						
U708-29	0010-20-42594 Accutor V, DPM 3, DPM 4, DPM 5, Passport V	Nellcor,non-Oximax	2.2M			B8
*U708-29M	M technology For old PM 6000 only	Nellcor,non-Oximax	2.2M			B8
*U708M-29	M technology For old PM 6000 only	Masimo LNCS	2.2M			B8
*U708-30M	M technology (0010-30-42625) DPM 3	Nellcor,non-Oximax	2.2M			B8
*U708M-30	DPM 4 DPM 5	Masimo LNCS	2.2M			B8
U708-22	0515-30-11221	Nellcor,non-Oximax	2.2M			C3
*U708-22M	M technology 0010-30-12452	non-Oximax Nellcor	2.2M			C3
*U708M-22	M technology 0010-30-12452	Masimo LNCS	2.2M			C3
U708-31	CSI Technology	CSI	2.2M			C3
*U708M-48	0010-30-42738 T5/T8 DPM 6, DPM 7	Masimo LNCS	2.2m			C3
U708-48	0010-20-42710 (C) T5/T8 DPM 6, DPM 7 Nellcor Non-Oximax	Nellcor,non-Oximax	2.2m			C3
U710X-48	0010-20-42712 T5/T8, DPM 6, DPM 7 Nellcor Oximax	Nellcor Oximax	3.0m			C4
*U708M-15	(0012-00-1652),1.2m (0012-00-1599),3.1m (0012-00-1653),4.3m Accutor V, Duo, Trio, Passport 2, Passport V, Spectrum/OR, V12, V21	Masimo LNCS				C3
U708-69	0012-00-1254 Trio, Passport 2, Spectrum	Nellcor,non-Oximax	2.2m			C3
U710-70	0012-00-1464 Accutor V, Duo, Trio, Passport 2, Passport V, Spectrum/OR, V12, V21	Nellcor Oximax	3.0m			C6
U710X-180	0010-20-42595 DPM4	Nellcor Oximax	3.0m			C3
Nellcor/Coviden 奈普特						
U708-69	SCP-10/MC-10	Nellcor, Non-Oximax	2.2m			C3
U710-85	M-200-13		3.0m			D7
U704-71	DEC-4		1.1m			B8
U708-71	DEC-8	Nellcor-Oximax	2.2m			B8
U710-70	DOC-10	Nellcor-Oximax	3.0m			C6
U710-70X	DOC-10 + Oximax adapter	Nellcor- Non Oximax	3.0m			C8

REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
Nellcor/Coviden 奈普特						
U701-01X U708-01X	Oximax adapter	Nellcor, Non-Oximax	0.3m 2.2m			C3
Nihon Kohden 光电						
U710-36	JL-900P BSM/PVM	NIHON KOHDEN	3.0m			C7
U710X-36	JL-650P	Oximax	3.0m			C7
U710M-36	JL-631P	M Tech	3.0m			C7
U710-106	JL-302T OPV-1500,OLV	NIHON KOHDEN	3.0m			C7
Nonin 燕牌						
U708-08		Nonin	2.2m			C3
Novamatrix 诺医						
U708-103	8853	Novamatrix	2.2m			C3
Pace tech/Generra						
U708-37		Nellcor non-Oximax	2.2m			B8
U708-38		Nellcor non-Oximax	2.2m			B8
U708-39		Nellcor non-Oximax	2.2m			C3
Philips/HP 飞利浦						
U708-40	M1940A/B	Philips	2.2m			C3
U708-41	M1941A/B		2.2m			C3
U708-20	M1900A/B	Nellcor non-Oximax	2.2m			C3
*U708M-20P	M1020-61102 CMS M1020A module.	Masimo LNOP	2.2m			C3
U708-43 U710-43	M1943A M1943AL	Nellcor non-Oximax	2.2m 3.0m			C3
U710X-43	M1943NL	Nellcor,Oximax	3.0m			C7
*U708-43M	M Technology	Nellcor non-Oximax	2.2m			C3
*U708M-43	989803148221 LNC MP10 (REF 2281)	Masimo LNCS sensor	2.2m			C3
*U708M-43P	M1020-61100	Masimo LNOP sensor	2.2m			C7

4 Accessories 配件

REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
Schiller席勒						
U708-17		BCI	2.2m			C3
*U708-17M	M tech.	Nellcor non-Oximax	2.2m			C3
*U708M-17	M tech.	Masimo LNCS sensor	2.2m			C3
U708-63 ODAM/Bruker		BCI	2.2m			C3
Siemens/Drager 西门子/德尔格						
U704-23	3368433	Nellcor non-Oximax	1.1m			C6
U708-23	3375834	Nellcor non-Oximax	2.2m			C6
U710X-23		Nellcor Oximax	3.0m			C7
*U704M-23P	MS18680	Masimo' LNOP	1.1m			C6
*U708M-23P	MS17041	Masimo' LNCS	2.2m			C6
*U704M-23	MS24303	Masimo' LNCS	1.1m			C6
*U708M-23	MS17522	Masimo' LNCS	2.2m			C6
U708-33	4533840	Nellcor non-Oximax	2.2m			C3
Spacelabs 太空						
U708-04	700-0002-00 90351-D/G, 90465, 90466, 90467, 90489,M77	Nellcor non-Oximax	2.2m			C3
U708-74	700-0030-00 90496,91496-U	Nellcor non-Oximax	2.2m			C6
U710X-74	700-0792-00 91496-N,91220	Nellcor-Oximax	3.0m			C7
*U708M-74	700-0906-00 91496-M	Masimo LNCS	2.2m			C6
*U708M-74P	700-0789-00 91496-M	Masimo LNOP	2.2m			C6
U701-04NV	175-0646-00 90351-D/G,90465, 90466,90467, 90489, 90651-08)	Novametrix- Hypertronic Con	0.3m			C3
U701-74NV	700-0029-00 (90496,91496-U)	Novametrix- Hypertronic Con	0.3m			C6
Welchallyn 伟伦						
U704X-104	001966-E (PIC 40TM Defibrillator PIC 50TM Defibrillator) 001965-U	Nellcor Oximax	2.2m			C7
U704-104	(PIC40 Defibrillator PIC50 Defibrillator)	Nellcor	1.1m			C7

*Refurbish Only



REF# 型号	Product Name 品名	Purpose of Use 用途	Package 包装	Price 价格	Figures 图标
Y04-E041	Ear sensor hanger 成人耳夹挂钩	Used with Unimed ear clip sensor U900 series or Other brands ear sensor(Gray) 用于尤迈900系列耳夹探头或其它品牌探头(灰)	1 pcs	A5	A
UEC-01	Adult ear clip Or as animal small size lingua clip 成人耳夹或作为小动物舌夹	Used with D-YS or Unimed U800 sensor(Black) 用与尤迈800系列或者D-YS探头(黑)	1 pcs	A5	B
UVC-01	Animal lingua clip,Big size 大号动物夹	Used with D-YS or Unimed U800 sensor(Gray) 用与尤迈800系列或者D-YS探头(灰)	1 pcs	A8	C
UVC-02	Animal lingua clip,Middle size 中号动物夹	Used with D-YS or Unimed U800 sensor(Gray) 用与尤迈800系列或者D-YS探头(灰)	1 pcs	A7	D
Y04-E019	Reusable Silicon Y sensor wrap tape 重复用Y型探头硅胶绑带	Used with Nellcor D-YS or Unimed U800 series sensors 用与尤迈800系列或者D-YS探头	1 pcs	A6	E
E023-12	Reusable 3 holes foam tape For All patient sizes 重复用3孔泡绵绑带	Used with Nellcor D-YS or Unimed U800 /U600 series sensors 用与尤迈800系列或者D-YS探头	12 pcs Per bag	B5/bag	F
E024-25	Disposable 3M adhesive textile tape For Adult and neonate 一次性3M无纺布绑带 (成人, 新生儿)	Used with Nellcor D-YS or Unimed U800 series sensors 用与尤迈800系列或者D-YS探头	25 pcs Per bag	B7/bag	G
E025-25	Disposable 3M adhesive textile tape for Pediatric and infant 一次性3M无纺布绑带 (儿童, 婴儿)	Used with Nellcor D-YS or Unimed U800 series sensors 用与尤迈800系列或者D-YS探头	25 pcs Per bag	B7/bag	H
T002-12	Velcro foam wrap tape 维克罗泡绵绑带	To fix the sensor cable onto the patient wrist 用于固定电缆在病人的手腕上	12 pcs per bag	C1.4/bag	I